



---

# User Manual

---



SCMS and CHAI would like to acknowledge the following organizations for their contributions and support:





---

# User Manual

---



# Contents

Chapter 1	
<b>Overview</b>	<b>1</b>
What ForLAB Can Do for You	1
Why Use ForLAB?	2
Who Should Use ForLAB?	3
Chapter 2	
<b>Getting Started</b>	<b>5</b>
System Requirements	5
How to Install ForLAB	5
Before You Begin	5
Installing ForLAB from a CD	6
How to Start ForLAB	6
Starting ForLAB from the Windows desktop	6
Starting ForLAB from the Windows taskbar	6
Reinstalling ForLAB	6
On-Line Help	7
Dashboards	7
Navigation Features	8
Tree view	8
Buttons	8
Scroll Bars	9
Pull-Down Menus	10
Keys for Editing Records	10
Chapter 3	
<b>Database Management</b>	<b>11</b>
Accessing the Database Control Panel	11
SQL Server Settings	11
Creating a New Database	12

<b>Switching between Databases</b>	<b>13</b>
<b>Backing up a Database</b>	<b>13</b>
<b>Restoring a Database</b>	<b>13</b>
<b>Chapter 4</b>	
<b>Test Profile</b>	<b>15</b>
<b>Testing Area vs. Testing Group vs. Test</b>	<b>15</b>
<b>Testing Area</b>	<b>15</b>
Adding Testing Area Data	16
Editing Testing Group Data	17
Testing Groups	17
TO ADD TESTING GROUP –	17
TO EDIT AN TESTING GROUP –	17
TO DELETE A TESTING GROUP –	17
Deleting a Testing Area	18
<b>Testing Groups</b>	<b>18</b>
Adding Testing Group Data	19
Editing Testing Group Data	19
Deleting a Testing Group	20
<b>Tests</b>	<b>20</b>
Adding Test Data	21
TO REMOVE A PRODUCT FROM A TEST—	21
Editing Test Data	22
Deleting a Test	22
<b>ART Tests</b>	<b>22</b>
Adding Test Data	23
TO REMOVE A REAGENT FROM A MORBIDIY TEST—	24
Adding Reagents for All Tests	25
TO REMOVE A REAGENT FROM AN ART TEST—	25
Editing Test Data	25
Deleting a Test	26

Chapter 5	
<b>Product Profile</b>	<b>27</b>
<b>Product Type</b>	<b>27</b>
Adding Product Type Data	27
Editing Product Type Data	28
Products	29
TO ADD PRODUCT –	29
TO EDIT AN PRODUCT –	29
TO DELETE A PRODUCT –	29
Deleting a Product Type	29
<b>Products</b>	<b>30</b>
Adding Products Data	30
Price	31
TO ADD A PRICE–	31
TO EDIT AN PRICE–	32
Editing Products Data	32
Deleting a Products	33
Chapter 6	
<b>Instruments</b>	<b>35</b>
<b>Instruments</b>	<b>35</b>
Adding Instrument Data	35
Editing Instrument Data	36
Deleting an Instrument	37
Chapter 7	
<b>Laboratory Profile</b>	<b>39</b>
<b>Regions/Districts/Provinces vs. Sites</b>	<b>39</b>
<b>Regions/Districts/Provinces</b>	<b>39</b>
Adding Region/District/Province Data	39
Editing Region/District/Province Data	40

<b>Sites</b>	<b>40</b>
TO ADD SITE –	40
TO EDIT AN SITE –	41
TO DELETE A SITE –	41
Deleting a Region/District/Province	41
<b>Site Category</b>	<b>41</b>
Adding Site Category Data	42
Editing Site Category Data	42
Deleting a Site Category	42
<b>Sites</b>	<b>43</b>
Adding Site Data	43
Referral Sites	44
TO ADD A REFERRAL SITE –	44
TO DELETE A REFERRAL SITE –	45
Site Status	45
TO ADD AN OPENING DATE –	45
TO EDIT AN OPENING DATE –	45
TO ADD A CLOSING DATE –	46
Instrument List	46
TO ADD AN INSTRUMENT–	46
TO DELETE AN INSTRUMENT–	47
Editing Site Data	47
Deleting a Site	48
<b>Chapter 8</b>	
<b>Morbidity Settings</b>	<b>49</b>
<b>Protocols</b>	<b>49</b>
Editing CD 4 Protocol Data	49
Editing Chemistry Protocol Data	50
Panels	51
TO ADD A PANEL –	51
TO EDIT A PANEL –	52
TO DELETE A PANEL –	52
Editing Hematology and Viral Load Protocol Data	53



<b>Panels</b>	<b>53</b>
TO ADD A PANEL –	53
TO EDIT A PANEL –	54
TO DELETE A PANEL –	54
<b>Editing Other Tests Protocol Data</b>	<b>54</b>
<b>Panels</b>	<b>55</b>
TO ADD A PANEL –	55
TO EDIT A PANEL –	56
TO DELETE A PANEL –	56
<b>Rapid Test Specification</b>	<b>56</b>
TO DEFINE THE SERIAL ALGORITHM—	58
TO DEFINE THE PARALLEL ALGORITHM—	58
<b>Chapter 9</b>	
<b>Forecasting: Consumption Methodology</b>	<b>61</b>
<b>Sources of Data</b>	<b>61</b>
<b>Consumption Methodology</b>	<b>61</b>
<b>Adding New Forecast</b>	<b>62</b>
ADD BY PRODUCT	64
ADD BY PRODUCTS AND SITES	65
IMPORT CONSUMPTION	66
<b>Editing Consumption Data</b>	<b>66</b>
TO ADD A REPORTING PERIOD TO THE BEGINNING OR END OF THE PERIODS LISTED—	67
TO REMOVE A REPORTING PERIOD FROM THE PERIODS LISTED –	67
<b>Calculating Forecast</b>	<b>67</b>
<b>Managing Forecast Result</b>	<b>69</b>
Forecast Summary	69
Printable Summary	70
<b>Chapter 10</b>	
<b>Forecasting: Service Statistics Methodology</b>	<b>71</b>
<b>Sources of Data</b>	<b>71</b>
<b>Service Statistics Methodology</b>	<b>71</b>
Adding New Forecast	72

ADD BY TEST	73
ADD BY TESTS AND SITES	74
IMPORT TESTING VOLUME	75
<b>Editing Service Statistics Data</b>	<b>75</b>
TO EDIT THE HISTORICAL TESTING VOLUME–	76
TO ADD A REPORTING PERIOD TO THE BEGINNING OR END OF THE PERIODS LISTED–	76
TO REMOVE A REPORTING PERIOD FROM THE PERIODS LISTED –	76
<b>Calculating Forecast</b>	<b>77</b>
<b>Managing Forecast Result</b>	<b>78</b>
Managing Forecast Result	78
Forecasted Service	79
Actual Service	79
Printable Summary	80
<b>Chapter 11</b>	
<b>Forecasting: Morbidity Methodology</b>	<b>83</b>
<b>Morbidity Methodology</b>	<b>83</b>
Adding New Forecast	83
Sites included in quantification	85
USING REGION AS A CATEGORY	86
USER DEFINED CATEGORY	86
Initial Patient Data	87
Current Patient Numbers by Site	87
OPTION 1	87
OPTION 2	88
Patient Treatment Target	89
Patient Ever Started on Treatment	89
EVER STARTED Patient Numbers by Site	90
OPTION 1	90
OPTION 2	91
OPTION 3	92
ART Patient Target	93
Site Target Calculator	94
OPTION 1	94

OPTION 2	94
OPTION 3	96
Site Specific Treatment Patient Numbers	97
Patient Pre-Treatment Targets	97
Site Pre-Treatment Target Calculator	99
OPTION 1	99
OPTION 2	99
OPTION 3	100
Site Specific Pre-Treatment Patient Numbers	101
Testing Efficiency	102
Testing Information	103
Adult Patient Behavior	104
Pediatric Patient Behavior	105
Rapid Test Algorithm by Site	106
Ordering and Inventory Assumptions	108
Reviewing Data	109
PATIENT BEHAVIOR	109
SERVICE AND PLATFORM	110
SUPPLY LIST	110
Editing Forecasts	111
Deleting Forecasts	111
<b>Calculating Results</b>	<b>112</b>
Patient Assumptions	112
CD4 Assumptions	112
Hematology, Viral Load, Chemistry, and Other Test Assumptions	112
HIV Rapid Test Assumptions	113
Demography Total Cost	113
Other Reports	113
TO NAVIGATE BETWEEN THE TWO CATEGORIES –	114
TO PRINT THE REPORTS –	114
TO VIEW THE PREVIOUS REPORTS –	114
Adjust Morbidity Supply Forecast	114



# Chapter 1

# Overview



*This section explains what ForLAB can do, why it should be used, and who can benefit from its use.*

The ForLAB Laboratory Quantification Tool (ForLAB), a software tool, was designed to help program managers forecast the needs for their laboratory products. ForLAB provides information needed to create and follow-up actions to ensure the regular and consistent stock of laboratory products at the program or national level. Consistency of stock is the first step in meeting the basic objective of any logistics system, which is to provide—

- the right quantities
- of the right commodities
- in the right condition
- in the right place
- at the right time
- for the right cost.

These are the six rights of logistics management.

## What ForLAB Can Do for You

ForLAB helps you achieve the right quantities at the right time.

ForLAB maintains—

- 1** Regions and facilities ordering tests including the number of working days for the laboratory and instruments at their disposal.
- 2** The test area, test group and test available for that instrument along with the product usage needed to perform the test.
- 3** Instruments available for each testing area along with the maximum throughput and number of control tests needed.
- 4** Products needed for each instrument along with the basic unit, minimum pack size and price for each.
- 5** Panels containing schedule information of patients that are recommended to receive tests, for each of the five major categories of testing (CD4, Chemistry, Hematology, Viral Load and other) along with other factors that may cause the actual volume of tests performed to differ significantly from the testing protocols like percentage of symptom directed test and test repeated due to clinical request.

With these data and an understanding of the methodologies required to calculate forecast consumption, ForLAB can calculate supply need based on the following methodologies—

- 1** Consumption methodology uses historic product consumption on the site level.
- 2** Service statistic methodology uses historic program service data at the site level.
- 3** Demographic methodology uses behavior of ART and pre-ART patients to forecast for a fixed 12 month period.

You can use this information with program policymakers, product suppliers, and donors to provide a rational basis for planning future product needs.

ForLAB is not the answer to every logistics question. It helps aggregate the quantity of each product needed to perform laboratory testing. ForLAB's utility is enhanced if your program has a standard operating procedure. Even without this underlying system, use ForLAB with whatever data are available. By beginning a rational and systematic forecasting process, you take the first step toward ensuring consistent stock levels.

## Why Use ForLAB?

Because of the large number of commodities that need to be forecasted, the process of forecasting the consumption of laboratory commodities requires multiple resources and can be time intensive. As a program manager, you face a complex forecasting environment, characterized by—

- 1** An extensive list of commodities dependent on the scope of the forecast and the level of standardization in the country doing the forecast.
- 2** Program targets may not consider the testing capacity at each facility or level of the system to provide services.
- 3** Lack of communication and coordination between policy makers, service providers, funding sources, and procurement agents on issues related to the selection, quantification, and procurement of laboratory supplies.
- 4** Multiple purposes of use in testing for most consumables making it challenging for staff to collect, aggregate, and report data that will be useful during quantification.

You need to review the data to determine its validity and usefulness for forecasting. Accuracy of the quantification will depend on the availability, completeness, and reliability of the data collected. To ensure credibility and ownership of the forecasting results, you should obtain consensus through a consultative and participatory process for gathering inputs and fostering discussion making among key stakeholders.

ForLAB can provide this information.

## Who Should Use ForLAB?

In a multi-product, multi-supplier environment, it is increasingly necessary that local program managers be empowered to do their own forecasting, pipeline monitoring, and procurement planning; they must also take charge of coordinating the activities of donors and local suppliers, as well as those of their own logistics management staff. Donor staff often has other priorities and little time to devote to the details of logistics management. Commercial suppliers have interests that may or may not correspond to the interests of your organization.

If you are the logistics manager or program manager for your organization, ForLAB can help.

While your managers and decision makers will be the primary users of ForLAB, the system can provide information to—

### **Purchasers/donors of commodities**

Staff who finance the purchase of commodities can use ForLAB reports and graphs to understand the future requirements of a commodity.

### **Host-country policymakers**

ForLAB reports and graphs can be used to help policymakers understand issues with the standard operating procedures and the implications of different decisions on the availability of the tests.





# Chapter 2

# Getting Started



*This section explains how to install ForLAB and presents some of its standard features.*

## System Requirements

The following resources are recommended for use with ForLAB—

CPU	Pentium III or higher
Operating System	Windows XP, Windows 2003 Server or higher
Memory	256 MB
Hard Drive Space	500 MB
Video Adapter	SVGA with 800x600 resolution
Miscellaneous	SQL Server 2005 or higher .NET framework 3.5 or higher Windows ReportView 2010 or higher

## How to Install ForLAB

ForLAB can be installed from a CD-ROM or the Internet.

### Before You Begin

Before installing ForLAB, you must install SQL Server 2005, .NET Framework 3.5, and Windows ReportView 2010 or higher. These files are provided with the ForLAB CD or can be downloaded from the following links:

**SQL Server 2005:** <http://www.microsoft.com/en-us/download/details.aspx?id=1695>

**.Net Framework 3.5:** <http://www.microsoft.com/en-us/download/details.aspx?id=22>

**Windows ReportView 2010:** <http://www.microsoft.com/en-us/download/details.aspx?id=6442>

## Installing ForLAB from a CD

- 1 Start Microsoft Windows.
- 2 Insert the ForLAB CD.

The ForLAB installation should begin automatically.

- 3 Follow the on-screen instructions.

If the installation does not begin automatically—

- 4 Click on Start on the Windows Taskbar.
- 5 Click on Run from the pop-up menu.
- 6 In the Command Line box, type x:setup (“x” is the letter of your CD-ROM drive).
- 7 Click on the OK button, and follow the on-screen instructions.

After ForLAB is successfully installed, the ForLAB shortcut will be displayed on your desktop.



## How to Start ForLAB

ForLAB can be started from the Windows desktop or the Windows taskbar.

### Starting ForLAB from the Windows desktop

From the Windows desktop—

- 1 Locate and double-click on the ForLAB icon to start the application.

### Starting ForLAB from the Windows taskbar

From the Windows taskbar—

- 1 Click on Start.
- 2 Click on Programs.
- 3 Locate and click on the ForLAB link.

## Reinstalling ForLAB

To reinstall ForLAB—

- 1 Place the ForLAB CD in your CD-ROM drive.
- 2 Start the install process, and follow the instructions on your screen.

During the process, a message box is displayed instructing you to remove ForLAB from your computer.

- 3 Click on the Remove button to remove ForLAB from your computer.
- 4 Click on the Finish button when prompted.

After ForLAB has been removed—

- 5 Repeat the ForLAB installation procedure.

## On-Line Help

ForLAB's On-Line Help is a supplement to the ForLAB User's Guide. On Line Help provides information about basic ForLAB functions.

- 1 Press the <F1> key to access On-Line Help.

## Dashboards

Throughout ForLAB you will find simple Dashboards (Figure 1) that visually display the information in the application. When first opening the applications, you will be taken to the main dashboard which displays six different graphs.

In addition to this dashboard, there are other graphs available when you click on each of the three profiles in the tree view, i.e. Test Profile, Product Profile, and Laboratory Profile. Each of the graphs in the tools has a feature to copy, save as image, print and zoom.

- 1 Right click on the chart to use these additional features.




Figure 1. Dashboard




Figure 2. Tree View

## Navigation Features

ForLAB has navigation features that help you move easily and quickly from one part of the system to another.

### Tree view

This feature offers a quick way to navigate through the system. The tree view is partially expanded when you start ForLAB, and it remains open on the left side of the screen. Figure 2 shows an example of the tree view.

To use the tree view—

- 1 Click on the option you need, to display its associated screen.

To collapse a group of tree view options—








- 2 Click on a tree view option to collapse the option and hide its sub-options.









To expand a group of tree view options—

- 3 Click on a tree view option to expand and show its sub-options.

## Buttons

The following buttons are most commonly used in ForLAB:

Button	Action
 Save and Close	Saves the record and closes the form.
 Save and New	Saves the record and opens a blank form.
 New...	Opens a blank form.
 Edit...	Opens the selected record.
 Delete	Deletes the selected record.
	Closes form without saving.
	Moves to first page.

Button	Action
	Moves to previous page.
	Moves to next page.
	Moves to last page.
	Edit record.
	Delete record.
	Adds current selection to the list.
	Saves the record and closes the form.
	Cancels changes and closes the form.

## Scroll Bars

Vertical and horizontal scroll bars are displayed on some ForLAB screens if there is more data than your screen can display at one time.

- 1 Click on one of the arrow buttons to move in the indicated direction:

Left	Up
Right	Down

Figure 3 displays an example of vertical and horizontal scroll bars.

**Tip:** To quickly scroll through a screen, click on the square button on the horizontal or vertical scroll bar, and hold down the left mouse button, while moving your mouse left and right or up and down.




Figure 3. Vertical and horizontal scroll bars

## Pull-Down Menus

**Tip:** You can also select an item from a pull-down menu by clicking on the menu and typing the first letter of the item you need.

Pull-down menus let you choose from a list. To use a pull-down menu—

- 1 Click on the arrow next to the menu to display its list.
- 2 Click on the item you want to select.

The selected item is displayed in the associated field.




Figure 4. Pull-down menu

## Keys for Editing Records

The following table lists keys you can use when editing records in ForLAB.

Keys for Editing Records	
F1	Displays On-Line Help.
Backspace	Deletes selection or character to the left of the cursor.
Delete	Deletes selection or character to the right of the cursor.
Esc	First press: Undoes changes to a field.
Insert	Second press: Undoes all changes to the record.
Tab	Toggles between insert and overwrite modes.
	Moves the cursor from one field to the next, and moves the cursor from one button to another without activating buttons.

## Chapter 3

# Database Management

*This section explains how to connect to and manage SQL Server databases in ForLAB.*

### Accessing the Database Control Panel

After you install ForLAB and click on the ForLAB icon for the first time, the Database Control Panel is displayed. (Figure 5) You can return to this screen by –

- 1 Click on the Data Management Menu Bar option.
- 2 Click on Database Control panel in the dropdown list.




Figure 5. Database Control Panel

### SQL Server Settings

After you install ForLAB and click on the ForLAB icon for the first time, the Database Control Panel is displayed. (Figure 5) You can return to this screen by –

- 1 Open the Database Control panel.
- 2 Click the Change button.

ForLAB displays the Connect to SQL Server screen similar to Figure 6.




Figure 6. Connect to SQL Server

- 3** Select the correct server from the Server dropdown list.
- 4** Select the Authentication method from the dropdown list.
- 5** Enter the appropriate User Name and Password for the authentication method selected.
- 6** Click Default to set this information as the default connection.
- 7** Click Quit to cancel changes and return to the Database Control Panel

Or, click Connect to connect to the selected server.

## Creating a New Database

To create a new database—

- 1** Open the Database Control Panel.
- 2** Click on the Create New Database button.
- 3** Type the database name you wish to call your new database.
- 4** Click OK.

A database with the name provided is created and displayed in the Databases list box.

- 5** Click Set as Default to make this database the default data file for ForLAB.
- 6** Click Quit to return to the application without changing databases or click Launch ForLAB to restart the application.

Each time you open ForLAB, the default database will automatically open.



## Switching between Databases

To switch to another database from within the ForLAB tool—

- 1 Open the Database Control panel.
- 2 Select the database you would like to open in the list box.
- 3 Click Set as Default to make this database the default data file for ForLAB.
- 4 Click Launch ForLAB to start the application.

## Backing up a Database

The Backup option of the ForLAB Database Control Panel allows you to make a copy of the current database. The copy can only be saved to your hard disk. Use this option to archive the data.

- 1 Open the Database Control panel.
- 2 Select the database you would like to open in the list box.
- 3 Click on the Backup button.

ForLAB displays a window so that you can select where the copy will be stored.

- 4 Select the location where you want to store the copy.
- 5 Click on the Open button to save the data to the selected location.

The data is saved.

## Restoring a Database

The Restore option of the ForLAB Database Control Panel allows you to make a restore a file not found in the list box. Use this option when sharing file with other users via email.

- 1 Open the Database Control panel.
- 2 Click on the Restore button.

ForLAB displays a window so that you can select the file to be restored.

- 3 Select the file you would like to restore.
- 4 Click on the Open button to add the database to the list box.
- 5 Click Set as Default to make this database the default data file for ForLAB.
- 6 Click Launch ForLAB to start the application.

The data is restored.



# Chapter 4

## Test Profile



*The Test Profile tree view options enable you to maintain the different Testing Areas, Testing Groups, and Tests available.*

### Testing Area vs. Testing Group vs. Test

A Testing Area is the platform of laboratory services that is to be incorporated in the quantification. For example cytology, serology, hematology, chemistry etc. Testing Group refers to the test profile that is associated with the testing group. For example for hematology, you have full blood count or in flow cytometry, you will have CD4 or in chemistry, you will have liver function tests, enzyme test etc.

### Testing Area

The Testing Area tree view option enables you to modify Testing Area names and other information.

- 1 Click on Settings tree view option.
- 2 Click on the General Settings tree view option.
- 3 Click on Test Profile tree view option.
- 4 Click on the Testing Area tree view option.

ForLAB opens the Testing Area screen (Figure 7).




Figure 7. Testing Area screen

## Adding Testing Area Data

With the Testing Area screen displayed—

- 1 Click on the Add New Testing Area link to display a blank Testing Area data entry form (Figure 8).


The screenshot shows a window titled "Testing Area" with a standard Windows-style title bar (minimize, maximize, close buttons). Inside the window, there are two buttons at the top: "Save and Close" (with a floppy disk icon) and "Save and New" (with a plus icon). Below these, there is a text input field for "Testing Area Name:". Underneath that is a checkbox labeled "Use In Demograph Methodology". Below the checkbox is a dropdown menu for "Class of Test:" which currently shows "CD4". Further down, there is a section titled "Testing Groups" in blue text. This section contains a table with a header "Group Name" and several empty rows. To the right of the table are three buttons: "New...", "Edit...", and "Delete".

Figure 8. Testing Area data entry form

- 2 Type the name of the Testing Area.

The Use in Demographic Methodology field specifies if this Testing Area can be included when forecasting based on Demographic Methodology.

- 3 Place a checkmark in the Use in Demographic Methodology field if Testing Area can be used.

The Class of Test field is only enabled when the Testing Area is included in Demographic Methodology.

- 4 Select the Class of Test from the dropdown field.

When you finish—

- 5 Click on the Save and Close button to save the data, and return to the Testing Area screen.

Or, click on the Save and New button to save the data, and display a blank Testing Area data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Testing Area screen.

## Editing Testing Group Data

Existing Testing Group data can be edited, as explained below.

With the Testing Group screen displayed—

- 1 Double click on the Testing Group in the Testing Group window that you want to edit.

Data associated with the selected Testing Area is displayed in the Testing Area data entry form.

- 2 Click on the field you want to edit, and enter the new data.

## Testing Groups

Testing Group Section lists the Testing Groups in the application that are included in the Testing Area. This section is only enabled when editing a Testing Area.

### TO ADD TESTING GROUP –

- 1 Click on the New button to display a blank Testing Group data entry form (Figure 10).
- 2 See Adding Testing Group Data on page 18 for information on completing the Testing Group data entry form.

### TO EDIT AN TESTING GROUP –

- 1 Select the Testing Group you want to edit from the list of Testing Groups.
- 2 Click the Edit button to display the Testing Group data entry form for that Testing Group.
- 3 See Adding Testing Group Data on page 18 for information on completing the Testing Group data entry form.

### TO DELETE A TESTING GROUP –

- 1 Select the Testing Group you want to delete from the list of Testing Groups.
- 2 Click the Delete button to display the Testing Group data entry form for that Testing Group.

If the Testing Group can be deleted, ForLAB will display a confirmation message.

- 3 Click on Yes to confirm the delete.

If you cannot delete the Testing Group because it is associated with other records ForLAB will display a message window.

When you finish—

- 3 Click on the Save and Close button to save the data, and return to the Testing Area screen.

Or, click on the Save and New button to save the data, and display a blank Testing Area data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Testing Area screen.

## Deleting a Testing Area

Use the Delete button to delete a Testing Area.

With Testing Area screen displayed—

- 1 Select the Testing Area that you want to delete.
- 2 Click on the Delete button.

If the Testing Area can be deleted, ForLAB will display a confirmation message.

- 3 Click on Yes to confirm the delete.

If you cannot delete the Testing Area because it is associated with other records ForLAB will display a message window.

## Testing Groups

The Testing Groups tree view option enables you to modify Testing Groups names and abbreviations.

- 1 Click on Settings tree view option.
- 2 Click on the General Settings tree view option.
- 3 Click on Test Profile tree view option.
- 4 Click on the Testing Groups tree view option.

ForLAB opens the Testing Groups screen (Figure 9).




Figure 9. Testing Group screen

## Adding Testing Group Data

With the Testing Groups screen displayed—

- 1 Click on the Add New Testing Group link to display a blank Testing Groups data entry form (Figure 10).

The figure shows a software window titled "Testing Group". At the top left are two buttons: "Save and Close" (with a floppy disk icon) and "Save and New" (with a plus icon). Below these are two input fields: "Group Name:" followed by a text box, and "Testing Area:" followed by a dropdown menu. Underneath is a section header "Tests" in blue. Below "Tests" is a table with one column labeled "Test Name". The table has several empty rows for data entry.

Figure 10. Testing Group data entry form

- 2 Type the name of the Testing Group.

The Testing Area field identifies the testing area where the Testing Groups belongs.

- 3 Click on the arrow next to the Testing Area field, and select an appropriate area from the pull-down menu.

The Tests field will list the entire test assigned to this Testing Group. When you finish—

- 4 Click on the Save and Close button to save the data, and return to the Testing Group screen.

Or, click on the Save and New button to save the data, and display a blank Testing Group data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Testing Group screen.

## Editing Testing Group Data

Existing Testing Group data can be edited, as explained below.

With the Testing Group screen displayed—

- 1 Double click on the Testing Group in the Testing Group window that you want to edit.

Data associated with the selected Testing Group is displayed in the Testing Group data entry form.

- 2 Click on the field you want to edit, and enter the new data.

When you finish—

- 3 Click on the Save and Close button to save the data, and return to the Testing Group screen.

Or, click on the Save and New button to save the data, and display a blank Testing Group data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Testing Group screen.

## Deleting a Testing Group

Use the Delete button to delete a Testing Group.

With the Testing Group screen displayed—

- 1 Select the Testing Group that you want to delete.
- 2 Click on the Delete button.

If the Testing Group can be deleted, ForLAB will display a confirmation message.

- 3 Click on Yes to confirm the delete.

If you cannot delete the Testing Group because it is associated with other records ForLAB will display a message window.

## Tests

The Tests tree view option enables you to modify Test names and abbreviations.

- 1 Click on Settings tree view option.
- 2 Click on the General Settings tree view option.
- 3 Click on Test Profile tree view option.
- 4 Click on the Tests tree view option.

ForLAB opens the Tests screen (Figure 11).




Figure 11. Tests screen



## Adding Test Data

With the Tests screen displayed—

- 1 Click on the Add New Test link to display a blank Test data entry form (Figure 12).

The image shows a software window titled "Test". At the top, there are two buttons: "Save and Close" and "Save and New". Below these are three input fields: "Test Name:" followed by a text box, "Testing Area:" followed by a dropdown menu, and "Testing Group:" followed by a dropdown menu. Below these fields is a section titled "Test Product Usage Rate". It contains two dropdown menus labeled "Ins." and "Pro.", followed by "Add" and "Remove" buttons. At the bottom of the window is a table with two columns: "Product" and "Usage Rate". The table is currently empty.

Figure 12. Test data entry form

- 2 Type the name of the Test.

The Testing Area field identifies the testing area where the Test belongs.

- 3 Click on the arrow next to the Testing Area field, and select an appropriate area from the pull-down menu.

The Testing Group field identifies the testing group where the Test belongs.

- 4 Click on the arrow next to the Testing Group field, and select an appropriate group from the pull-down menu.

The Product Usage field allows the user to select the instrument and product used to perform this Test along with the Usage rate for the product.

- 5 Click on the arrow next to the Ins. field, and select an appropriate instrument from the pull-down menu.
- 6 Click on the arrow next to the Pro. field, and select an appropriate Product from the pull-down menu.
- 7 Click on Add button to add the Instrument/Product to the Test.
- 8 Double Click the Usage Rate field to enter the correct usage rate for the product.

### TO REMOVE A PRODUCT FROM A TEST—

- 1 Select the Product from the Product Usage list.
- 2 Click on Remove button to remove the Product from the Test.

When you finish—

- 9 Click on the Save and Close button to save the data, and return to the Test screen.

Or, click on the Save and New button to save the data, and display a blank Test data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Test screen.

## Editing Test Data

Existing Test data can be edited, as explained below.

With the Test screen displayed—

- 1 Double click on the Test in the Test window that you want to edit.

Data associated with the selected Test is displayed in the Test data entry form.

- 2 Click on the field you want to edit, and enter the new data.

When you finish—

- 3 Click on the Save and Close button to save the data, and return to the Test screen.

Or, click on the Save and New button to save the data, and display a blank Test data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Test screen.

## Deleting a Test

Use the Delete button to delete a Test.

With the Test screen displayed—

- 1 Select the Test that you want to delete.
- 2 Click on the Delete button.

If the Test can be deleted, ForLAB will display a confirmation message.

- 3 Click on Yes to confirm the delete.

If you cannot delete the Test because it is associated with other records ForLAB will display a message window.

## ART Tests

ART Tests are necessary since forecasting with the morbidity methodology models behavior of ART and pre-ART patients.

- 1 Click on Settings tree view option.
- 2 Click on the General Settings tree view option.
- 3 Click on Test Profile tree view option.

- 4 Click on the ART Tests tree view option to expand the list of ART Test types.
- 5 Click on the ART test type

ForLAB opens the ART Tests screen for the selected ART test type (Figure 13).




Figure 13. ART Tests screen

## Adding Test Data

With the ART Tests screen displayed—

- Click on the Add New Test link to display a blank Morbidity Test data entry form (Figure 14).




Figure 14. Morbidity Test data entry form

The Test Name field is populated when you save the test.

The Platform field identifies the instrument used in the test.

- 1 Click on the arrow next to the Platform field, and select an appropriate instrument from the pull-down menu.
- 2 Click Save and Close.

The Selected Reagent According to Test field allows the user to select the reagents used to perform this Test along with the Usage rate for the product and the process for calculating the quantification (Quantify According to...).

- 3 Double click on the test you just added.
- 4 Click Add Reagent to add a new reagent to Selected Reagent According to Test (Figure 15).

Figure 15. Select Reagents

- 5 Click on the arrow next to the Quantify According to... field, and select the appropriate process from the pull-down menu.
- 6 Click on the Product/Reagent Name that is will be quantified in this test
- 7 Click the Select button.

The selected Product/Reagent Names will be added to the list at the bottom of the form.

- 8 Repeat steps 6 -8 for each Reagent quantified in this test.
- 9 When finished, click OK to save the form and return to the Morbidity Test data entry form.

Or, click on the Cancel button to cancel all changes and return to the Morbidity Test data entry form.

#### TO REMOVE A REAGENT FROM A MORBIDIY TEST—

- 1 Select the Product from the Selected Reagent According to Test list.
  - 2 Click on Remove Reagent link to remove the Product from the Test.
- When you finish—
- 10 Click on the Save and Close button to save the data, and return to the Test screen.

Or, click on the Save and New button to save the data, and display a blank Morbidity Test data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the ART Test screen.

## Adding Reagents for All Tests

In some cases, products can be used in an ART test independent of the specific instrument used. These products can be entered separately. With the ART Tests screen displayed—

- 1 Click Select Product Link to add a new reagent to the ART Test (Figure 15).
- 2 Click on the arrow next to the Quantify According to... field, and select the appropriate process from the pull-down menu.
- 3 Click on the Product/Reagent Name that is will be quantified in this test.
- 4 Click the Select button.

The selected Product/Reagent Names will be added to the list at the bottom of the form.

- 5 Repeat steps 2 -4 for additional Reagents.
- 6 When finished, click OK to save the form and return to the ART Test screen.

Or, click on the Cancel button to cancel all changes and return to the ART Test screen.

### TO REMOVE A REAGENT FROM AN ART TEST—

- 1 Select the Product from the list.
- 2 Click on Remove link to remove the Product from the ART Test screen.

## Editing Test Data

Existing ART Test data can be edited, as explained below.

With the appropriate ART Test screen displayed—

- 1 Double click on the Test in the Test window that you want to edit.

Data associated with the selected Test is displayed in the Morbidity Test data entry form.

- 2 Click on the field you want to edit, and enter the new data.

When you finish—

- 3 Click on the Save and Close button to save the data, and return to the ART Test screen.

Or, click on the Save and New button to save the data, and display a blank Morbidity Test data entry form

Or, click on the red X in the upper right hand corner to disregard all changes and return to the ART Test screen.

## Deleting a Test

Use the Delete button to delete a Test.

With the ART Test screen displayed—

- 1** Select the Test that you want to delete.
- 2** Click on the Delete button.

If the Test can be deleted, ForLAB will display a confirmation message.

- 3** Click on Yes to confirm the delete.

If you cannot delete the Test because it is associated with other records ForLAB will display a message window.

# Chapter 5

## Product Profile



*The Product Profile tree view options enable you to maintain the different Product Types and Products available.*

### Product Type

The Product Type tree view option enables you to modify Product Type names and other information.

- 1 Click on Settings tree view option.
- 2 Click on the General Settings tree view option.
- 3 Click on Product Profile tree view option.
- 4 Click on the Product Type tree view option.

ForLAB opens the Product Type screen (Figure 16).




Figure 16. Product Type screen

### Adding Product Type Data

With the Product Type screen displayed—

- 1 Click on the Add New Product Type link to display a blank Product Type data entry form (Figure 17).

Figure 17. Product Type data entry screen

- 2 Type the name of the Product Type.
- 3 Type a description of the Product Type.

The Use in Demographic Methodology field specifies if this Product Type can be included when forecasting based on Demographic Methodology.

- 4 Place a checkmark in the Use in Demographic Methodology field if Product Type can be used.

The Class of Test field is only enabled when the Product Type is included in Demographic Methodology.

- 5 Select the Class of Test from the dropdown field.

When you finish—

- 6 Click on the Save and Close button to save the data, and return to the Product Type screen (Figure 16).

Or, click on the Save and New button to save the data, and display a blank Product Type screen.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Product Type screen.

## Editing Product Type Data

Existing Product Type data can be edited, as explained below.

With the Product Type screen displayed—

- 1 Double click on the Product Type in the Product Type window that you want to edit.

Data associated with the selected Product Type is displayed in the Product Type data entry form.

- 2 Click on the field you want to edit, and enter the new data.



## Products

Product Section lists the Products in the application that are included in the Product Type. This section is only enabled when editing a Product Type.

### TO ADD PRODUCT –

- 1 Click on the New button to display a blank Product data entry form (Figure 19).
- 2 See Adding Product Data on page 18 for information on completing the Product data entry form.

### TO EDIT AN PRODUCT –

- 1 Select the Product you want to edit from the list of Product.
- 2 Click the Edit button to display the Product data entry form for that Testing Group.
- 3 See Adding Product Data on page 18 for information on completing the Product data entry form.

### TO DELETE A PRODUCT –

- 1 Select the Product you want to delete from the list of Product.
- 2 Click the Delete button to display the Product data entry form for that Product.

If the Product can be deleted, ForLAB will display a confirmation message.

- 3 Click on Yes to confirm the delete.

If you cannot delete the Product because it is associated with other records ForLAB will display a message window.

When you finish—

- 3 Click on the Save and Close button to save the data, and return to the Product Type screen.

Or, click on the Save and New button to save the data, and display a blank Product Type data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Product Type screen.

## Deleting a Product Type

Use the Delete button to delete a Product Type.

With Product Type screen displayed—

- 1 Select the Product Type that you want to delete.
- 2 Click on the Delete button.

If the Product Type can be deleted, ForLAB will display a confirmation message.

- 3 Click on Yes to confirm the delete.

If you cannot delete the Product Type because it is associated with other records ForLAB will display a message window.

## Products

The Products tree view option enables you to modify Products names and abbreviations.

- 1 Click on Settings tree view option.
- 2 Click on the General Settings tree view option.
- 3 Click on Product Profile tree view option.
- 4 Click on the Products tree view option.

ForLAB opens the Testing Groups screen (Figure 18).



Figure 18. Product Screen

## Adding Products Data

With the Products screen displayed—

- 1 Click on the Add New Product link to display a blank Product data entry form (Figure 19).

Price	Pack Size	As of Date

Figure 19. Product data entry form

- 2 Type the name of the Products.

The Product Type field identifies the product type where the Products belong.

- 3 Click on the arrow next to the Product Type field, and select an appropriate area from the pull-down menu.

The Serial Number field identifies the serial number of the product used in procurement.

- 4 Type the serial number of the Product.

The Specification field identifies the manufacturer's specification or the national specification of the instrument (i.e. Blood Collection Tube, 4 ml, plastic with screw cap, red top, and 100/pack). Note that this specification is important in identifying the products during procurement. If this specification is wrong, you may have the wrong product procured for you. It is therefore in the interest of the country that the specification is complete. If you know the manufacturers specification, you can use it.

- 5 Type the specifications of the Product.

The Basic Unit is the unit of measure used to count the product at the usage level.

- 6 Type the basic unit of the Product.

The Minimum Packs/Site is the minimum pack size available at the central store that is used for procurement or issued to the lower level, if the same.

- 7 Type the minimum packs per site for the Product.

The Notes field identifies optional notes on the product.

- 8 Type optional notes for the Product.

The Test Specification field indicates which category the test belongs if it is part of broader tests. For example ALT is part of enzyme tests. But this can be left blank for tests that are not part of broader specification. For example CD4 is a known test.


- 9 Type the test specification for the Product.

## Price

The Price group identifies the price, pack size and effective date for the product.

### TO ADD A PRICE--

- 1 Click on New button to open the Price Input dialog box (Figure 20).



The image shows a 'Price input' dialog box with three input fields: 'Price', 'Pack Size', and 'As of Date'. The 'As of Date' field is pre-filled with '10/23/2013' and has a calendar icon to its right. At the bottom of the dialog are 'Ok' and 'Cancel' buttons.

Figure 20. Price Input dialog box

- 2** Type in the Price of the Product.
- 3** Type in the Pack size of the Product.
- 4** Type in the As Of Date of the Price.
- 5** Click on the OK button to save the data, and return to the Product data entry screen.

Or, click on the Cancel button to return to the Product data entry screen without saving.

#### **TO EDIT AN PRICE–**

- 1** Select the price that is no longer available at the site
- 2** Click on the Edit Button to open the Price Input dialog box (Figure 20).
- 3** Click on the field you want to edit, and enter the new data.
- 4** Click on the OK button to save the data, and return to the Product data entry screen.

Or, click on the Cancel button to return to the Product data entry screen without saving.

When you finish—

- 10** Click on the Save and Close button to save the data, and return to the Products screen.

Or, click on the Save and New button to save the data, and display a blank Product data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Product screen.

## **Editing Products Data**

Existing Products data can be edited, as explained below.

With the Products screen displayed—

- 1** Double click on the Product in the Products window that you want to edit.

Data associated with the selected Product is displayed in the Product data entry form.

- 2** Click on the field you want to edit, and enter the new data.

When you finish—

- 3** Click on the Save and Close button to save the data, and return to the Products screen.

Or, click on the Save and New button to save the data, and display a blank Product data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Product screen.

## Deleting a Products

Use the Delete button to delete a Product.

With the Products screen displayed—

- 1** Select the Product that you want to delete.
- 2** Click on the Delete button.

If the Product can be deleted, ForLAB will display a confirmation message.

- 3** Click on Yes to confirm the delete.

If you cannot delete the Product because it is associated with other records ForLAB will display a message window.



# Chapter 6

# Instruments



*The Instruments tree view options enable you to maintain the different Instruments available.*

## Instruments

The Instruments tree view option enables you to modify Instrument names and other information.

- 1 Click on Settings tree view option.
- 2 Click on the General Settings tree view option.
- 3 Click on the Instruments tree view option.

ForLAB opens the Instruments screen (Figure 21).




Figure 21. Instrument screen

## Adding Instrument Data

With the Instruments screen displayed—

- 1 Click on the Add New Instrument link to display a blank Instrument data entry form (Figure 22).

Figure 22. Instrument data entry form

- 2 Type the name of the Instrument.

The Testing Area field identifies the testing area where the Instrument belongs.

- 3 Click on the arrow next to the Testing Area field, and select an appropriate area from the pull-down menu.

The Max Throughput field identifies the maximum number of test an instrument can process per hour, day, or month. It is important to note that the max throughput is correctly identify if you are using test per hour or test per month to determine days or hours work to calculate utilization. This is not the same as the number of specimens, as many tests can be done on a single specimen.

- 4 Type in the maximum throughput for the instrument.

The Control Test Needed field identifies the number of control tests needed to be run and how often.

- 5 Type in the number of control test to be run for each applicable time period.

When you finish—

- 6 Click on the Save and Close button to save the data, and return to the Instruments screen.

Or, click on the Save and New button to save the data, and display a blank Instrument data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Instrument screen.

## Editing Instrument Data

Existing Instrument data can be edited, as explained below.

With the Instrument screen displayed—

- 1 Double click on the Instrument in the Instrument window that you want to edit.



Data associated with the selected Instrument is displayed in the Instrument data entry form.

- 2 Click on the field you want to edit, and enter the new data.

## **Deleting an Instrument**

Use the Delete button to delete an Instrument.

With Instrument screen displayed—

- 1 Select the Instrument that you want to delete.
- 2 Click on the Delete button.

If the Instrument can be deleted, ForLAB will display a confirmation message.


- 3 Click on Yes to confirm the delete.

If you cannot delete the Instrument because it is associated with other records ForLAB will display a message window.



# Chapter 7

# Laboratory Profile



*The Laboratory Profile tree view options enable you to maintain the different Regions/Districts/Provinces within the country along with the sites within those regions.*

## Regions/Districts/Provinces vs. Sites

A Region/District/Province is an administrative division or unit of a country. Within each Region/District/Province there are multiple facilities providing laboratory services. Sites refer to the individual facility providing and referring patients for laboratory services.

## Regions/Districts/Provinces

The Regions/Districts/Provinces tree view option enables you to modify Region/District/Province names and abbreviations along with assigning sites within the region.

- 1 Click on the Regions/Districts/Provinces tree view option.

ForLAB opens the Regions/Districts/Provinces screen (Figure 23).




Figure 23. Regions/Districts/Provinces screen

## Adding Region/District/Province Data

With the Region/District/Province screen displayed—

- 1 Click on the Add New Region/District/Province link to display a blank Region/District/Province data entry form (Figure 24).

Figure 24. Region / District / Province data entry form

- 2 Type the name of the Region/District Province.

The Short Name field holds the Region/District/Provinces short name, which is used to help identify the Region/District/Province.

- 3 Type the short name of the Region/District/Province.

When you finish—

- 4 Click on the Save and Close button to save the data, and return to the Region/District/Province screen.

Or, click on the Save and New button to save the data, and display a blank Region/District/Province data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Region/District/Province screen.

## Editing Region/District/Province Data

Existing Region/District/Province data can be edited, as explained below.

With the Region/District/Province screen displayed—

- 1 Double click on the Region/District/Province in the Region/District/Province window that you want to edit.

Data associated with the selected Region/District/Province is displayed in the Region/District/Province data entry form.

- 2 Click on the field you want to edit, and enter the new data.

## Sites

Site Section lists the sites in the application that reside in the Region/District/Province. This section is only enabled when editing a Region/District/Province.

### TO ADD SITE –

- 1 Click on the New button to display a blank Site data entry form (Figure 28).
- 2 See Adding Site Data on page 18 for information on completing the Site data entry form.

**TO EDIT AN SITE –**

- 1** Select the Site you want to edit from the list of sites.
- 2** Click the Edit button to display the Site data entry form for that site.
- 3** See Adding Site Data on page 18 for information on completing the Site data entry form.

**TO DELETE A SITE –**

- 1** Select the Site you want to delete from the list of sites.
- 2** Click the Delete button to display the Site data entry form for that site.
- 3** If the Site can be deleted, ForLAB will display a confirmation message.
- 4** Click on Yes to confirm the delete.

If you cannot delete the Site because it is associated with other records ForLAB will display a message window.

When you finish—

- 3** Click on the Save and Close button to save the data, and return to the Region/District/Province screen.

Or, click on the Save and New button to save the data, and display a blank Region/District/Province data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Region/District/Province screen.

## **Deleting a Region/District/Province**

Use the Delete button to delete a Region/District/Province.

With the Region/District/Province screen displayed—

- 1** Select the Region/District/Province that you want to delete.
- 2** Click on the Delete button.

If the Region/District/Province can be deleted, ForLAB will display a confirmation message.

- 3** Click on Yes to confirm the delete.

If you cannot delete the Region/District/Province because it is associated with other records ForLAB will display a message window.

## **Site Category**

The Site Categories tree view option enables you to modify Site Category names and abbreviations.

- 1** Click on the Site Categories tree view option.

ForLAB opens the Site Categories screen (Figure 25).




Figure 25. Site Category screen

## Adding Site Category Data

With the Site Category screen displayed—

- 1 Click on the Add New Site Category link to display a blank Site Category data entry form (Figure 26).




Figure 26. Site Category data entry form

- 2 Type the name of the Site Category.

When you finish—

- 3 Click on the OK button to save the data, and return to the Site Category screen.

Or, click on the Cancel button to return to the Site Category screen without saving.

## Editing Site Category Data

Existing Site Category data can be edited, as explained below.

With the Site Category screen displayed—

- 1 Double click on the Site Category in the Site Category window that you want to edit.

Data associated with the selected Site Category is displayed in the Site Category data entry form.

- 2 Enter the new data.

When you finish—

- 3 Click on the Save button to save the data, and return to the Site Category screen.

## Deleting a Site Category

Use the Delete button to delete a Site Category.

With the Site Category screen displayed—

- 1 Select the Site Category that you want to delete.
- 2 Click on the Delete button.

If the Site Category can be deleted, ForLAB will display a confirmation message.

- 3 Click on Yes to confirm the delete.

If you cannot delete the Site Category because it is associated with other records ForLAB will display a message window.

## Sites

The Sites tree view option enables you to modify Site names and abbreviations.

- 1 Click on the Sites tree view option.

ForLAB opens the Sites screen (Figure 27).




Figure 27. Sites screen

## Adding Site Data

With the Site screen displayed—

- 1 Click on the Add New Site link to display a blank Site data entry form (Figure 28).

Figure 28. Site Data entry form

**2** Type the name of the Site.

The Region/District/Province field identifies the region, district or province where the site resides.

**3** Click on the arrow next to the Region/District/Province field, and select an appropriate area from the pull-down menu.

The Site Category field identifies the primary group to which the new site belongs.

**4** Click on the arrow next to the Site Category field, and select an appropriate category from the pull-down menu.

The General Working Days\Month field identifies the number of days per month that the site is operating. The default value for this field is 22 days.

**5** Type in the number of workings days per month.

The Site Level field identifies the facility level to which the new site belongs.

**6** Click on the arrow next to the Site Level field, and select an appropriate level from the pull-down menu.

The Test Days per Month section identifies the number of testing days in each of the five diagnostic areas (CD4, Chemistry, Hematology, Viral Load, and Other). Note that the value entered for any one of these diagnostic areas cannot be greater than the number of working days per month entered previously.

**7** Type the number of testing days for diagnostic area at the new site.

## Referral Sites

The Referral Sites section identifies the site tests are referred to for the diagnostic area if the new site is not providing the service.

### TO ADD A REFERRAL SITE –

- 1** Click the ellipse button for the diagnostic area whose service is not provided by the site. (Figure 29).



- 2 Click on the arrow next to the Region/District/Province field, and select an appropriate area from the pull-down menu.




Figure 29. Add Referral Site

- 3 Click on the arrow next to the Referral Site field, and select an appropriate Site from the pull-down menu.

**Tip:** Only sites that have instruments for the diagnostic area will be included in the pull-down menu

When you finish—

- 4 Click Add button to add the Referral Site and return to the Site data entry form.

Or, click Cancel to cancel these entries and return to the Site data entry form.

#### TO DELETE A REFERRAL SITE –

- 1 Click on Del button for the appropriate diagnostic area.

## Site Status

Site Status indicates if the new site is currently operational.

#### TO ADD AN OPENING DATE –

- 1 Click on Open link to display the Opening Date dialog box (Figure 30).




Figure 30. Opening Date dialog box

- 2 Click on the arrow next to the Opening Date field, and select an appropriate date from the calendar.

When you finish—

- 3 Click on the OK button to save the data, and return to the Site screen

Or, click Red X to cancel these entries and return to the Site data entry form.

#### TO EDIT AN OPENING DATE –

- 1 Click on Edit link to display the Opening Date dialog box.

Data associated with the selected Site Status is displayed in the Opening Data dialog box.

- 2 Click on the arrow next to the Opening Date field, and select an appropriate date from the calendar.

When you finish—

- 3 Click on the OK button to save the data, and return to the Site screen  
Or, click Red X to cancel these entries and return to the Site data entry form.

#### TO ADD A CLOSING DATE –

- 1 Click on Close link to display the Closing Date dialog box (Figure 31).




Figure 31. Closing Date dialog box

- 2 Click on the arrow next to the Closing Date field, and select an appropriate date from the calendar.

When you finish—

- 3 Click on the OK button to save the data, and return to the Site screen  
Or, click Red X to cancel these entries and return to the Site data entry form.

## Instrument List

The Instrument List shows the complete list of diagnostic instruments available at the new site.

#### TO ADD AN INSTRUMENT–

- 1 Click on Add Instrument link to open the Select Instrument dialog box (Figure 32).




Figure 32. Select Instrument dialog box

- 2 Select each instrument available at the site. (To select more than one instrument, hold the Ctrl key while selecting instruments.)
- 3 Click Select to add the instrument(s) to the site and return to the Site data entry screen.

To enter the Quantity of the instrument at the site –

- 4** Click in the Quantify column for the instrument you want to change.
- 5** Type in the Quantity of that instrument at the site
- 6** Hit Enter

% Test Run is the percentage of times that instrument is used in the testing area. The sum of % Test Run for instruments under the same testing area should be 100%. To enter the % Test Run—

- 7** Click in the % Test Run column for the instrument you want to change.
- 8** Type in the percentage of Tests Run on that instrument at the site
- 9** Hit Enter

#### **TO DELETE AN INSTRUMENT—**

- 1** Select the instrument that is no longer available at the site
- 2** Click on the Remove link.

When you finish—

- 3** Click on the Save and Close button to save the data, and return to the Site screen.

Or, click on the Save and New button to save the data, and display a blank Site data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Site screen.

## **Editing Site Data**

Existing Site data can be edited, as explained below.

With the Site screen displayed—

- 1** Double click on the Site in the Site window that you want to edit.

Data associated with the selected Site is displayed in the Site data entry form.

- 2** Click on the field you want to edit, and enter the new data.

When you finish—

- 3** Click on the Save and Close button to save the data, and return to the Site screen.

Or, click on the Save and New button to save the data, and display a blank Site data entry form.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Site screen.

## Deleting a Site

Use the Delete button to delete a Site.

With the Site screen displayed—

- 1** Select the Site that you want to delete.
- 2** Click on the Delete button.

If the Site can be deleted, ForLAB will display a confirmation message.

- 3** Click on Yes to confirm the delete.

If you cannot delete the Site because it is associated with other records ForLAB will display a message window.

# Chapter 8

# Morbidity Settings

*The Morbidity Setting tree view options enable you to easily prepare for forecasting using the Demographic/Morbidity methodology.*

## Protocols

The Protocols tree view options enable you to maintain the different schedules that patients are recommended to receive tests. Understanding the testing protocol is critical in estimating how many tests a patient population will receive, and therefore how many testing supplies will be required.

The Protocol tree view option enables you to modify Protocols.

- 1 Click on Settings tree view option.
- 2 Click on Morbidity Settings tree view option.
- 3 Click on the Protocols tree view option.

ForLAB opens the Protocols screen (Figure 33). The four major categories of testing protocols as well as a fifth category for other protocols are displayed.




Figure 33. Protocols screen

## Editing CD 4 Protocol Data

With the Protocols screen displayed—

- 1 Double click on the CD4 Protocol to display the CD4 Protocol Settings form (Figure 34).

	Month1	Month2	Month3	Month4	Month5	Month6	Month7	Month8	Month9	Month10	Tests/Year after month 12
Adults in Treatment	0	0	0	0	0	0	0	0	0	0	
Pediatrics in Treatment	0	0	0	0	0	0	0	0	0	0	
Adults pre-ART	0	0	0	0	0	0	0	0	0	0	
Pediatrics pre-ART	0	0	0	0	0	0	0	0	0	0	

Figure 34. CD4 Protocol Setting

- 2 Type the Description of the protocol.

The % of test repeated field indicates the percentage of tests that are repeated.

- 3 Type in the % of test repeated for the protocol.

The % symptom directed test/year field indicates the percentage of patients each year who receive a symptom directed test, in addition to the testing schedule entered below.

- 4 Type in the % symptom-directed test/year for the protocol.

The Patient protocol as number of blood samples drawn for tests in a given month (after the initial test) fields indicate the number of times that each patient group should receive this panel each month during the first year of initiating treatment.

- 5 Type in the number of blood samples drawn for each patient group for each month (1 -12.)

The Tests/year after month 12 indicates the number of times that each patient group should receive this panel each year after the first year of treatment.

- 6 Type in the number of blood samples drawn for each patient group for each year after the first year of treatment.

When you finish—

- 7 Click on the Save and Close button to save the data, and return to the Protocols screen.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Protocols screen.

## Editing Chemistry Protocol Data

With the Protocols screen displayed—

- 1 Double click on the Chemistry Protocol to display the Chemistry Protocol Settings form (Figure 35).

Figure 35. Chemistry Protocol Setting

**2** Type the Description of the protocol.

The % of test repeated field indicates the percentage of tests that are repeated.

**3** Type in the % of test repeated for the protocol.

## Panels

The Panel tab lists the different panels applied.

### TO ADD A PANEL –

**1** Click on the New button to display a blank Protocol Panel data entry form (Figure 36).

Figure 36. Chemistry Protocol Panel

The Panel Name field describes what circumstances this panel will apply.

**2** Type in the Panel Name.

The Test to Include field indicates which tests are included in the panel.

**3** Place a check mark in the box for each test included in the panel.

The Patient protocol as number of blood samples drawn for tests in a given month fields indicate the number of times that each patient group should receive this panel each month.

**4** Type in the % of New Patients on the Panel for each patient group.

**5** Type in the number of blood samples drawn for each patient group for each month (1 -12.)

- 6** Type in the % of Pre-Existing Patient on the Panel for each patient group.
- 7** Type in the number of blood samples drawn for each patient group for each year after the first year of treatment.

#### TO EDIT A PANEL –

- 1** Select the Panel you want to edit from the list of Panels.
- 2** Click the Edit button to display the Protocol Panel data entry form for that Panel.
- 3** Click on the field you want to edit, and enter the new data.
- 4** Click on the OK button to save the data, and return to the Protocol Settings screen.

Or, click on the Cancel button to return to the Protocol Settings screen without saving.

#### TO DELETE A PANEL –

- 1** Select the Panel you want to delete from the list of Panels.
- 2** Click the Delete button to display the Protocol Panel data entry form for that Panel.

If the Panel can be deleted, ForLAB will display a confirmation message.

- 3** Click on Yes to confirm the delete.

If you cannot delete the Panel because it is associated with other records ForLAB will display a message window.

Panel	% Total patients requiring 1 symptom-directed test/year	ALT	AST	CHO	GLC	CRE	TG	GGT	ALP	AMY	CO2
Adults in Treatment	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Pediatrics in Treatment	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Adults pre-ART	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0
Pediatrics pre-ART	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0	0.0

Figure 37. Chemistry Protocol Setting

The % symptom directed test/year tab (Figure 37) indicates the percentage of patients each year who receive a symptom directed test, in addition to the testing schedule entered for each patient group.

- 4** Type in the % symptom-directed test/year for each patient group and test in the panel.

When you finish—

- 5** Click on the Save and Close button to save the data, and return to the Protocols screen.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Protocols screen.



## Editing Hematology and Viral Load Protocol Data

With the Protocols screen displayed—

- 1 Double click on the Hematology or Viral Load Protocol to display the Hematology or Viral Load Protocol Settings form (Figure 38).

Figure 38. Hematology/Viral Load Protocol Setting

- 2 Type the Description of the protocol.

The % of test repeated field indicates the percentage of tests that are repeated.

- 3 Type in the % of test repeated for the protocol.

The % symptom directed test/year field indicates the percentage of patients each year who receive a symptom directed test, in addition to the testing schedule entered below.

- 4 Type in the % symptom-directed test/year for the protocol.

## Panels

The Panel tab lists the different panels applied.

### TO ADD A PANEL –

- 1 Click on the New button to display a blank Protocol Panel data entry form (Figure 39).

Figure 39. Hematology/Viral Load Protocol Panel

The Panel Name field describes what circumstances this panel will apply.

- 2 Type in the Panel Name.

The Patient protocol as number of blood samples drawn for tests in a given

month fields indicate the number of times that each patient group should receive this panel each month.

- 3** Type in the % of New Patients on the Panel for each patient group.
- 4** Type in the number of blood samples drawn for each patient group for each month (1 -12.)
- 5** Type in the % of Pre-Existing Patient on the Panel for each patient group.
- 6** Type in the number of blood samples drawn for each patient group for each year after the first year of treatment.

#### **TO EDIT A PANEL –**

- 1** Select the Panel you want to edit from the list of Panels.
- 2** Click the Edit button to display the Protocol Panel data entry form for that Panel.
- 3** Click on the field you want to edit, and enter the new data.
- 4** Click on the OK button to save the data, and return to the Protocol Settings screen.

Or, click on the Cancel button to return to the Protocol Settings screen without saving.

#### **TO DELETE A PANEL –**

- 1** Select the Panel you want to delete from the list of Panels.
- 2** Click the Delete button to display the Protocol Panel data entry form for that Panel.

If the Panel can be deleted, ForLAB will display a confirmation message.

- 3** Click on Yes to confirm the delete.

If you cannot delete the Panel because it is associated with other records ForLAB will display a message window.

When you finish—

- 4** Click on the Save and Close button to save the data, and return to the Protocols screen.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Protocols screen.

## **Editing Other Tests Protocol Data**

With the Protocols screen displayed—

- 1** Double click on the Other Tests Protocol to display the Other Tests Protocol Settings form (Figure 40).




Figure 40. Other Tests Protocol Settings

- 2 Type the Description of the protocol.
- 3 The % of test repeated field indicates the percentage of tests that are repeated.
- 4 Type in the % of test repeated for the protocol.

Panels

The Panel tab lists the different panels applied.

TO ADD A PANEL –

- 1 Click on the New button to display a blank Protocol Panel data entry form (Error! Reference source not found.).




Figure 41. Other Tests Protocol Panel

The Panel Name field describes what circumstances this panel will apply.

- 2 Type in the Panel Name.

The Test to Include field indicates which tests are included in the panel.

- 3 Place a check mark in the box for each test included in the panel.

The Patient protocol as number of blood samples drawn for tests in a given month fields indicate the number of times that each patient group should receive this panel each month.

- 4 Type in the % of New Patients on the Panel for each patient group.
- 5 Type in the number of blood samples drawn for each patient group for

each month (1 -12.)

- 6** Type in the % of Pre-Existing Patient on the Panel for each patient group.
- 7** Type in the number of blood samples drawn for each patient group for each year after the first year of treatment.

#### **TO EDIT A PANEL –**

- 1** Select the Panel you want to edit from the list of Panels.
- 2** Click the Edit button to display the Protocol Panel data entry form for that Panel.
- 3** Click on the field you want to edit, and enter the new data.
- 4** Click on the OK button to save the data, and return to the Protocol Settings screen.

Or, click on the Cancel button to return to the Protocol Settings screen without saving.

#### **TO DELETE A PANEL –**

- 1** Select the Panel you want to delete from the list of Panels.
- 2** Click the Delete button to display the Protocol Panel data entry form for that Panel.

If the Panel can be deleted, ForLAB will display a confirmation message.

- 3** Click on Yes to confirm the delete.

If you cannot delete the Panel because it is associated with other records ForLAB will display a message window.

The % symptom directed test/year tab indicates the percentage of patients each year who receive a symptom directed test, in addition to the testing schedule entered for each patient group.

- 5** Type in the % symptom-directed test/year for each patient group and test in the panel.

When you finish—

- 6** Click on the Save and Close button to save the data, and return to the Protocols screen.

Or, click on the red X in the upper right hand corner to disregard all changes and return to the Protocols screen.

## **Rapid Test Specification**

Before specifying the rapid test settings for each site in the quantification the rapid test specifications should be defined for each rapid test in the system.




Figure 42. Protocols screen

Rapid test specification includes sensitivity and specificity. Sensitivity refers to the percentage of actual HIV that the test detects, and specificity refers to the percentage of positive detections that are actually HIV. In other words, sensitivity measures the likelihood that a negative result is a true negative, and specificity measures the likelihood that a positive result is a true positive.

- 1 Click on Settings tree view option.
- 2 Click on Morbidity Settings tree view option.
- 3 Click on Rapid Test Specification.

The Rapid Test Specifications form will be displayed.(Figure 43)

HIV rapid test	Test Sensitivity	Test Specificity	HIV + False Negative	HIV - False Positive
<b>Screening</b>				
Determine HIV 1&2 with or without buffer	99.8	99.8	0.1	0.1
	0	0	100	100
	0	0	100	100
<b>Confirmatory</b>				
Test	99.7	99.8	0.3	0.2
	0	0	100	100
	0	0	100	100
<b>Tie-Breaker</b>				
Stat Pak HIV 1/2, 20 Tests	99.7	99.9	0.3	0.1
	0	0	100	100
	0	0	100	100

Figure 43. Rapid Test Specification screen

- 4 Using the drop down list, select up to three Screening Tests.
- Similarly, select up to three Confirmatory Tests.

And, select up to three Tie-Breaker Tests.

- 5 Click Save Rapid Test button to populate the Serial and Parallel Tabs with the tests selected.

#### TO DEFINE THE SERIAL ALGORITHM—

- 1 Click on Serial tab.
- 2 Double click in the Test Sensitivity cell.
- 3 Type in the value for the Test Sensitivity
- 4 Double click in the Test Specificity cell.
- 5 Type in the value for the Test Specificity.
- 6 Repeat steps 2 to 5 for each test.

#### TO DEFINE THE PARALLEL ALGORITHM—

- 1 Click on the Parallel tab.

**Rapid Test Specifications**

**Screening Tests**

Rapid Product

Determine VIH kit/100 with chase buffer

---NONE---

---NONE---

**Confirmatory Tests**

Rapid Product

Test

---NONE---

---NONE---

**Tie-Breaker Tests**

Tie-Breaker Product

Stat-Pak HIV 1/2, 20 Tests

---NONE---

---NONE---

Save Rapid Test

**Serial** **Parallel**

HIV rapid test	Test Sensitivity	Test Specificity	HIV + False Negative	HIV - False Positive
<b>Screening</b>				
Determine VIH kit/100 with chase buffer	0	0	100	100
	0	0	100	100
	0	0	100	100
<b>Confirmatory</b>				
Test	0	0	100	100
	0	0	100	100
	0	0	100	100
<b>Tie-Breaker</b>				
Stat-Pak HIV 1/2, 20 Tests	0	0	100	100
	0	0	100	100
	0	0	100	100

**Rules for Progressing Through Algorithm**

Both Negative

Stop

Both Positive

Stop

Discordant

Proceed

Figure 44. Rapid Test Specification, parallel tab

- 2 Double click in the Test Sensitivity cell.
- 3 Type in the value for the Test Sensitivity
- 4 Double click in the Test Specificity cell.
- 5 Type in the value for the Test Specificity.
- 6 Repeat steps 2 to 5 for each test.

The Rules for Progressing through Algorithm define the action to take for each of the possible results of the first two parallel tests, i.e. both positive, both negative, and discordant. For each result, there are only two possible actions, to stop or to proceed with another test.

- 7** Select the appropriate action for Both Positive.
  - 8** Select the appropriate action for Both Negative.
  - 9** Select the appropriate action for Discordant.
- 5** Click on the Save Rapid Test button.





## Chapter 9

# Forecasting: Consumption Methodology



*The Consumption method uses historical data on the use of laboratory commodities to calculate the quantities that will be demanded in the future.*

### Sources of Data

Consumption data are historical data on the actual quantities of a product that have been dispensed to patients or used at a service delivery point within a given time period, and are typically reported per month or per quarter. Daily consumption data can be found in pharmacy dispensing registers, laboratory registers, or other point of service registers. Where a well-functioning LMIS captures and aggregates these data from service delivery points, aggregated consumption data can be found in monthly and annual facility-level and program-level reports. For antiretroviral drugs (ARVs), consumption data would be the actual quantity of each ARV dispensed to ART patients. For HIV tests, consumption data or “usage data” are the actual number of HIV tests used over a given period. For laboratory supplies, consumption data are the actual number of laboratory commodities used.

### Consumption Methodology

The Consumption Methodology tree view option enables you to modify forecast information.

- 1 Click on Quantification Process tree view option.
- 2 Click on the Consumption Methodology tree view option.

ForLAB opens the Consumption Forecast Methodology screen (Figure 45) that will show the list of forecasts that have been entered into the database.




Figure 45. Consumption Forecast Methodology screen

## Adding New Forecast

With the Consumption Forecast Methodology screen displayed—

- 1 Click on the Add New Forecast link to display a blank Quantification Process data entry form (Figure 46).

Figure 46. Quantification Process data entry form

The Forecast ID field specifies the unique id and short text describing the forecast.

- 2** Type the Forecast ID.

The Scope of the Forecast field categorizes the forecasts to help the user to have different categories of forecasts. This field can be Custom, National, or Global.

- 3** Select the Scope of the Forecast from the drop down list.

The Reporting Period field is independent on the type of historical data collected for the forecast. It can be bimonthly, monthly, quarterly, or yearly.

- 4** Select the Reporting Period from the dropdown field.

The Start Date field is the date the forecast period will start.

- 5** Type in the Start Date for the forecast.

The Forecasting Period field is the number of periods we wanted to include in the forecast.

- 6** Type in the Forecasting Period.

The Forecast Date field and Last Modified Date field are system dates on which the forecast was defined and modified. Neither of these fields is editable.

The Consumption Data Usage field indicates the way service data is to be entered and used in the system. Options include:

- Consumption Data available for all sites – use this option if you have data reported from each site in the forecast
- Consumption Data available for reported sites – use this option if you only have a portion of sites reported and you want to then use this data for non-reported sites.
- Consumption Data available for site category – use this option if you want to categorize sites in groups based on level and testing behaviors and enter service data for the group.

- 3** Click in the radial button for the appropriate data usage.

When you finish—

- 4** Click on the Save Forecast Info button to save the data.

- 5** Enter Consumption Data

With the Quantification Process screen displayed—

- 6** Click on the Consumption tab in the Quantification Process window that you want to edit (Figure 47).

Figure 47. Quantification Process data entry form

There are three ways to add Service Statistic Data.

#### ADD BY PRODUCT

The List of Sites field lists all the sites in the system. The Selected Sites field lists all the sites included in the forecast defined.

- 1 Select the sites in the List of Sites field to be included in the forecast.
- 2 Click on Select Site link to add to add selected sites to the Selected Sites list.

Or click Select All link to add all sites to the Selected Sites list.

To remove a Selected Site from the Selected Sites list—

- 1 Select the sites in the Selected Sites field to be removed from the forecast.
- 2 Click on Remove Site link to add to remove the sites from the Selected Sites list.

The Add by Product link allows users to add a product for the selected site with the number of reporting periods provided. To add the products—

- 3 Select the Site in the Selected Sites list.
- 4 Click on the Add By Product link.

The Select Product form is opened (Figure 48). This form is used to add a product for the selected site with the number of reporting periods provided.

Product Name	Serial No	Basic Unit	Specification

Figure 48. Select Product

- 5** Select the products to add to the forecast. (Select more than one product by holding down the ctrl button while clicking on the product.)
  - 6** Type in the number of Reporting Periods for the products.
  - 7** Click Select button to add the selected products to the selected sites.
- Or, click cancel to return to the Quantification Process screen.

#### ADD BY PRODUCTS AND SITES

The Add Test and Sites link allows the user to add Product and sites simultaneously with the number of reporting periods provided. To add the Products and sites—

- 1** Click on the Add by Products and Sites link.

The Select Product form is opened (Figure 49). This form is used to add a product for the selected sites with the number of reporting periods provided.

Product Name	Serial No	Basic Unit	Specification
BD FACS Count-Control Kit	SN002	test	
BD FACS Count-FacClean	SN003	litre	
BD FACS Count-FacFlow	SN005	litre	
BD FACS Count-FacRinse	SN004	litre	
BD FACS Count-Thermal Pa	SN008	ml	
Determine VIH - kit/100 wt.	SN167	test	
FACS Count CD3/4 Reage	SN001	test	
Fully (bio)system Agulle	SN038	Agulle	
Fully (bio)system Calibrator S	SN034	ml	
Fully (bio)system Cleaner Re	SN040	ml	
Fully (bio)system Control Ser	SN032	ml	
Fully (bio)system Control Ser	SN033	ml	
Fully (bio)system Creatinine	SN028	ml	
Fully (bio)system Cuvettes, B	SN035	cuvette	
Fully (bio)system Glucose 1x	SN029	ml	
Fully (bio)system Lampe	SN037	lampe	
Fully (bio)system Paper impri	SN039	rouleau	
Fully (bio)system Sample cu	SN036	cup	
Fully (bio)system TGO/ASA	SN030	ml	
Fully (bio)system TGP/ALT	SN031	ml	
SD Boline HIV 1/2 3.0	SN165	test	
Symex Ki 21 - Cell clean C	SN158	ml	
Symex Ki 21 - CellPack P	SN156	litre	
Symex Ki 21 - Control Eigh	SN159	Tube	
Symex Ki 21 - Control Fig	SN161	Tube	

Figure 49. Select Products form

- 2** Select the sites to add to the forecast by placing a check in the box next to each site.
- 3** Select the products to add to the forecast. (Select more than one product by holding down the ctrl button while clicking on the product.)

- 4 Type in the number of Reporting Periods for the products.
  - 5 Click Select button to add the selected products to the selected sites.
- Or, click cancel to return to the Quantification Process screen.

### IMPORT CONSUMPTION

The Import Consumption link allows users to import bulk consumption from the excel template. To import the tests and sites—

- 1 Click on the Import Consumption link.

The Import form is opened (Figure 50).

Figure 50. Import form

- 2 Click the Browse button to show the Open File dialog box.
- 3 Find the excel template file and click open.

The file path is populated on the import form.

- 4 Click Import.
  - 5 Double click Consumption sheet in the list.
  - 6 Click Save to save the data to the Quantification form
- Or click clear to clear out the data and select another file to import.

### Editing Consumption Data

Once a test is added to the site using one of the methods above, the Product Consumption will populate and can be updated directly (Figure 51).

Figure 51. Product Consumption data entry form

- 1 Select a site in the Selected Site list.
- 2 Click on the cell you want to update and enter the correct value.

#### TO ADD A REPORTING PERIOD TO THE BEGINNING OR END OF THE PERIODS LISTED–

- 1 Click Add Actual Consumption link.

The popup window displays the start date of the forecast.

- 2 Enter the number of reporting periods to add.
- 3 Click OK.

#### TO REMOVE A REPORTING PERIOD FROM THE PERIODS LISTED –

- 1 Click on the reporting period to be deleted.

ForLAB will only allow period at the beginning or end of the forecast to be delete.

- 2 Click Remove Consumption link.

When you finish—

- 3 Click on the Save Consumption button to save the data.
- 4 Click on the red X to close the form and return to the Consumption Forecast Methodology screen.

**Note:** ForLAB will not allow the user to delete service if only 3 records

## Calculating Forecast

With the Consumption Forecast Methodology screen displayed—

- 1 Click on Forecast Button at the top left corner of the application.

The Forecasting Tools form will display. (Figure 52) The Filter Group allows user to filter forecasts in the system using methodology and data usage.

Figure 52. Forecasting Tools

- 2 Select Consumption Methodology from dropdown list.
- 3 Select Data Usage from the dropdown list.

Once the ID of the Forecast is selected “Do Forecast” button will be enabled. Under Forecasting Parameters users can change the regression of the forecast and the Wastage and Add by % to be added on the forecast result.

- 4 Select the Regression Type from the dropdown list.
- 5 Type in the Wastage %.
- 6 Type in the Program Growth % over the forecast period.
- 7 Click on Do Forecast button to start the Forecast.

The Forecasting process takes time to finalize; when complete the Forecast Result Form is shown. (Figure 53)

Figure 53. Forecast Result

This form displays the forecasting parameters entered on the previous form along with any messages. Forecast accuracy is measured using the Mean Absolute Percent Error (MAPE) and is defined as the average of percentage errors. The MAPE Chart tab (Figure 55) will display this indicator as a chart and the MAPE Result tab (Figure 54) will display the calculated percentage errors. The results are also color coded to quickly determine under forecasts, over forecasts and insufficient data in the forecast.






Figure 55. MAPE Chart tab

Product Name	June 2012	September 2012	December 2012	April 2013
FACS Count CD3	7.00	11.00	35.00	26.00
BD FACS Count...	1.00	0.17	-1.33	-2.67
BD Test CD3/...	5.83	-8.67	-5.57	8.33
BD Calibre 3 Be...	0.00	0.00	0.00	0.00
BD FACS Lysing...	2.00	-1.67	-3.50	-4.60
Easy CD4/CD4%	4.00	-0.67	-5.33	-10.67
Fully Biosystem	6.00	-1.00	-8.00	-16.00
Cobas C111 Part...	0.00	0.00	0.00	0.00

Figure 56. MAPE Result Tab

## Managing Forecast Result

Once a forecast is done, different summary information of the result can be seen on the forecast form. (Figure 56) From this form, results can be viewed for each site included in the forecast or as all sites combined.

## Forecast Summary




Figure 56. Forecast Summary

This report shows the forecasted cost summary and cost trend for the overall forecast or for a specific selected site.

- 1 Select all or an individual Site in the Select Site to view Chart field.

The Forecast Summary will refresh based on the selection.

## Printable Summary

To view detailed printable version of a forecast result –




Figure 57. *frmReportViewer*

- 1 Click on Show Summary button.

The frmReportViewer form will display. (Figure 57)

- 2 Click on printer button to send report to printer

Or, Click on disk button to send report to PDF, Word, or Excel.

## Chapter 10

# Forecasting: Service Statistics Methodology

*The Service Statistics method uses historical numbers of tests performed over a given period of time to estimate future testing demands.*

## Sources of Data


Service statistics data are historical, program-level or facility-level data on the number of patient visits to facilities, the number of services provided, or the number of people who received a specific service or treatment within a given time period. Service statistics data can be found in program monitoring reports, HMIS data, facility-level data on service utilization and attendance rates, or in patient records. In some programs, the LMIS captures a limited number of service statistics. For ARVs, service statistics data would be the total number of ART patients on treatment at a facility, or perhaps the total number of patient visits to a facility at a point in time. For HIV tests, service statistics would be the total number of clients tested during a certain period. For laboratory supplies, service statistics are the total number of tests performed in a certain period (e.g., CD4 count tests performed in a given quarter).

## Service Statistics Methodology

The Service Statistics Methodology tree view option enables you to modify forecast information.

- 1 Click on Quantification Process tree view option.
- 2 Click on the Service Statistics Methodology tree view option.

ForLAB opens the Service Statistics Forecast Methodology screen (Figure 58) that will show the list of forecasts that have been entered into the database.



Status	Forecast ID	Scope	Forecast Date	Stat Date	Reporting Period	Forecast Period	Method	Weightage	Last Updated
CLD	service/f	NATO	10/29/2013	10/1/2	Monthly	3	Linear	0.00	10/29/2013

Figure 58. Service Statistics Forecast Methodology screen

## Adding New Forecast

With the Service Statistics Forecast Methodology screen displayed—

- 1 Click on the Add New Forecast link to display a blank Quantification Process data entry form (Figure 59).

Figure 59. Quantification Process data entry form

The Forecast ID field specifies the unique id and short text describing the forecast.

- 2 Type the Forecast ID.

The Scope of the Forecast field categorizes the forecasts to help the user to have different categories of forecasts. This field can be Custom, National, or Global.

- 3 Select the Scope of the Forecast from the drop down list.

The Reporting Period field is independent on the type of historical data collected for the forecast. It can be bimonthly, monthly, quarterly, or yearly.

- 4 Select the Reporting Period from the dropdown field.

The Start Date field is the date the forecast period will start.

- 5 Type in the Start Date for the forecast.

The Forecasting Period field is the number of periods we wanted to include in the forecast.

- 6 Type in the Forecasting Period.

The Forecast Date field and Last Modified Date field are system dates on which the forecast was defined and modified. Neither of these fields is editable.

The Service Statistics Data Usage field indicates the way service data is to be entered and used in the system. Options include:

- When you finish—

- ## Enter Service Statistics Data

- Click on the Service Statistics tab in the Quantification Process window that you want to edit (Figure 60).

### ADD BY TEST

- 1 Select the sites in the List of Sites field to be included in the forecast.



- Or click **Select All** link to add all sites to the **Selected Sites** list.

- 1 Select the sites in the Selected Sites field to be removed from the forecast.

- The Add by Test link allows users to add a test for the selected site with the number of reporting periods provided. To add the tests—

**3** Select the Site in the Selected Sites list.

**4** Click on the Add by Test link.

The Select Test form is opened (Figure 61). This form is used to add a test for the selected site with the number of reporting periods provided.

Figure 61. Select Test Form

**5** Select the tests to add to the forecast. (Select more than one test by holding down the ctrl button while clicking on the test.)

**6** Type in the number of Reporting Periods for the tests.

**7** Click Select button to add the selected tests to the selected sites.

Or, click cancel to return to the Quantification Process screen.

#### ADD BY TESTS AND SITES

The Add Test and Sites link allows the user to add test and sites simultaneously with the number of reporting periods provided. To add the tests and sites—

**1** Click on the Add by Tests and Sites link.

The Select Test form is opened (Figure 62). This form is used to add a test for the selected sites with the number of reporting periods provided.

Figure 62. Select Test Form

**2** Select the sites to add to the forecast by placing a check in the box next to each site.

- 3** Select the tests to add to the forecast. (Select more than one test by holding down the ctrl button while clicking on the test.)
  - 4** Type in the number of Reporting Periods for the tests.
  - 5** Click Select button to add the selected tests to the selected sites.
- Or, click cancel to return to the Quantification Process screen.

### IMPORT TESTING VOLUME

The Import Testing Volume link allows users to import bulk service data from the excel template. To import the tests and sites—

- 1** Click on the Import Testing Volume link.

The Import form is opened (Figure 63).

Figure 63. *Import form*

- 2** Click the Browse button to show the Open File dialog box.
- 3** Find the excel template file and click open.

The file path is populated on the import form.

- 4** Click Import.
  - 5** Double click Service Statistics sheet in the list.
  - 6** Click Save to save the data to the Quantification form
- Or click clear to clear out the data and select another file to import.

## Editing Service Statistics Data

Once a test is added to the site using one of the methods above, the Historical Testing Volume will populate and can be updated directly (Figure 64).



**Quantification Process**

General Service Statistics

**List of Sites** [Select All](#) [Select Site](#)

Site Name

**Selected Sites** [Remove Site](#)

Edoras Lab  
Mona Medical Center  
Ongliath Labs

**Historical Testing Volume** [Import Testing Volume](#) [Add Test and Sites](#) [Add By Test](#) [Add Test V](#) [Remove Test V](#)

May-2013 June-20... July-2013 August... Septem... October...

**Glucose**

Amount	45.00	53.00	62.00	72.00	83.00	97.00
StockOut	1	0	0	0	0	0
Instrum	1	0	0	0	0	0
Adjusted	49.50	53.00	62.00	72.00	83.00	97.00
Note	Adjusted	-	-	-	-	-

**ALT**

Amount	46.00	52.00	61.00	70.00	82.00	95.00
StockOut	0	0	0	0	0	0
Instrum	0	0	0	0	0	0
Adjusted	46.00	52.00	61.00	70.00	82.00	95.00
Note	-	-	-	-	-	-

**Hgb**

Amount	50.00	58.00	67.00	78.00	90.00	104.00
StockOut	1	0	0	0	0	0
Instrum	1	0	0	0	0	0
Adjusted	55.00	58.00	67.00	78.00	90.00	104.00
Note	Adjusted	-	-	-	-	-

**CD4**

Amount	53.00	61.00	71.00	83.00	96.00	111.00
StockOut	1	0	0	0	0	0
Instrum	1	0	0	0	0	0
Adjusted	58.30	61.00	71.00	83.00	96.00	111.00
Note	Adjusted	-	-	-	-	-

**Creatinine**

Active Site:  Total Sites:

Test Count:  Total Test Count:

[Save Consumption](#)

Figure 64. Historical Testing Volume data entry form

**TO EDIT THE HISTORICAL TESTING VOLUME–**

- 1 Select a site in the Selected Site list.
- 2 Click on the cell you want to update and enter the correct value.

**TO ADD A REPORTING PERIOD TO THE BEGINNING OR END OF THE PERIODS LISTED–**

- 1 Click Add Test V link.

The popup window displays the start date of the forecast.

- 2 Enter the number of reporting periods to add.
- 3 Click OK.

**TO REMOVE A REPORTING PERIOD FROM THE PERIODS LISTED –**

- 1 Click on the reporting period to be deleted.

ForLAB will only allow period at the beginning or end of the forecast to be delete.

- 2 Click Remove Test V link.

When you finish—

- 3 Click on the Save Consumption button to save the data.
- 4 Click on the red X to close the form and return to the Service Statistics Forecast Methodology screen.

**Note:** ForLAB will not allow the user to delete service if only 3 records



## Calculating Forecast

With the Service Statistics Forecast Methodology screen displayed—

- 1 Click on Forecast Button at the top left corner of the application.

The Forecasting Tools form will display. (Figure 65) The Filter Group allows user to filter forecasts in the system using methodology and data usage.

Figure 65. Forecasting Tools

- 2 Select Service Statistics Methodology from dropdown list.
- 3 Select Data Usage from the dropdown list.

Once the ID of the Forecast is selected “Do Forecast” button will be enabled. Under Forecasting Parameters users can change the regression of the forecast and the Wastage and Add by % to be added on the forecast result.

- 4 Select the Regression Type from the dropdown list.
- 5 Type in the Wastage %.
- 6 Type in the Program Growth % over the forecast period.
- 7 Click on Do Forecast button to start the Forecast.

The Forecasting process takes time to finalize; when complete the Forecast Result Form is shown. (Figure 67)

Figure 67. Forecast Result

This form displays the forecasting parameters entered on the previous form along with any messages. Forecast accuracy is measured using the Mean Absolute Percent Error (MAPE) and is defined as the average of percentage errors. The MAPE Chart tab (Figure 68) will display this indicator as a chart and the MAPE

Result tab (Figure 66) will display the calculated percentage errors. The results are also color coded to quickly determine under forecasts, over forecasts and insufficient data in the forecast.




Figure 68. MAPE Chart tab



Figure 69. MAPE Result Tab

## Managing Forecast Result

Once a forecast is done, different summary information of the result can be seen on the forecast form. (Figure 69) From this form, results can be viewed for each site included in the forecast or as all sites combined.

### Managing Forecast Result

Once a forecast is done, different summary information of the result can be seen on the forecast form. (Figure 69) From this form, results can be viewed for each site included in the forecast or as all sites combined.




Figure 69. Forecast Summary

- 1 Select all or an individual Site in the Select Site to view Chart field.

The Forecast Summary will refresh based on the selection

## Forecasted Service

This report shows the forecasted diagnostics contribution, Instrument utilization and instrument distribution for the instruments under the overall forecast or for a specific selected site.




Figure 70. Forecasted Service

## Actual Service

This report shows the actual diagnostics contribution, Instrument utilization and instrument distribution for the instruments under the overall forecast or for a specific selected site.




Figure 71--Actual Service

## Printable Summary

To view detailed printable version of a forecast result –

- 1 Click on Show Summary button.

The frmReportViewer form will display. (Figure 72)




Figure 72. frmReportViewer

2 Click on printer button to send report to printer

Or, Click on disk button to send report to PDF, Word, or Excel.



## Chapter 11

# Forecasting: Morbidity Methodology



*The morbidity method forecasts the quantity of laboratory commodities demanded for the diagnosis, care and treatment of specific diseases based on existing disease prevalence and program targets.*

## Morbidity Methodology

The Morbidity Methodology tree view option enables you to modify forecast information.

- 1 Click on Quantification Process tree view option.
- 2 Click on the Morbidity Methodology tree view option.

ForLAB opens the Morbidity Forecast Methodology screen (Figure 73) that will show the list of forecasts that have been entered into the database.




Figure 73. Morbidity Forecast Methodology screen

## Adding New Forecast

With the Morbidity Forecast Methodology screen displayed—

- 1 Click on the Add New Forecast link to display a blank Morbidity Quantification data entry form (Figure 74).

The Title of Quantification field specifies a short title describing the forecast.

Figure 74. Morbidity Quantification

- 2 Type the Title of Quantification.
- 3 Select the Date of Quantification.

The Description, Purpose and Key Assumptions field specifies the description of the quantification, including the purpose and key assumptions that have been made to create the quantification.

- 4 Type the Description, Purpose and Key Assumptions.

The Forecast Start Date field is the date the forecast period will start.

- 5 Type in the Forecast Start Date for the forecast.

It is important to understand the distinction between the forecast timeframe and the budget period. The forecast timeframe always includes the twelve months starting with the forecast start date. However, it is sometimes desirable to perform quantification for only a specific number of months within that 12-month forecast timeframe. This period is the Budget Period. The Budget Start field is disabled since this will always be the same as the Forecast Start month. The Budget End field is the last month in to include in the budget. ForLAB will include all months from the Budget Start through the Budget End.

- 6 Select the Budget End.

The Type of Algorithm field specifies the desired algorithm for rapid testing to be used in the forecast. Serial indicates that each patient is tested once initially, those that receive a positive result are retested to confirm the result, and patients with discordant results are tested a third time. Parallel indicates that all patients are tested twice initially and patients with discordant results have the option to be tested a third time.

- 7 Click in the radial button for the appropriate algorithm.

When you finish—



- 8 Click on the Save button to save the data.

Upon saving the data, the working area to the right will be activated (Figure 75). At the top of this working area are three buttons and a description of the Working Area page. These buttons will move you through the working area. The Next button will move you to the next page in the working area. The Back button will return you to the previous page. The Go Home button will return you to the first page in the working area. Please remember to use the Save button in the top left hand corner to save your work as you go through the working area.

Figure 75. *Morbidity Quantification with Working Area enabled*

## Sites included in quantification

The home screen in the working area is the Sites Included in Quantification screen. (Figure 76)

Figure 76. *Sites included in quantification*

The Categorize Site By field allows the user to categorize sites when adding them into the forecast. There are two ways of categorizing sites: (1) using region as a category; or (2) defining your own custom category.

- 1 Click the Use Region as Categories or Use User Defined Categories radial button.

#### USING REGION AS A CATEGORY

- 1 Click in the Use region as categories radial button.

The complete list of regions will display in the List of Categories field.

#### USER DEFINED CATEGORY

- 1 Click in the Use user defined categories radial button.

Repeat the following steps 2 through 4 for each user defined category needed:

- 2 Click Add New button.

The Category Name field describes the category being defined.

- 3 Edit Category Name.

- 4 Click Save to save the user defined category

Or, click Delete to remove it from the list.

The complete list of categories will display in the List of Categories field.

- 2 Select a category in the List of Categories.

- 3 Click Add Site link to add sites to the quantification for the specified category.

The Select Site form will display (Figure 77).

Region/District/Province Name	Site Name	Type
DS ABENGOUROU	BETHESDA ABENGOUROU	
DS ABENGOUROU	CAT ABENGOUROU	
DS ABENGOUROU	CHR ABENGOUROU	
DS ABENGOUROU	CSU AFFALIKRO	
DS ABENGOUROU	CSU AMELEKIA	
DS ABENGOUROU	CSU APPROMPRON	
DS ABENGOUROU	CSU BETTIE	
DS ABENGOUROU	CSU NIABLE	
DS ABENGOUROU	CSU SANKADIOKRO	
DS ABENGOUROU	CSU YAKASSE-FEYASSE	
DS ABENGOUROU	CSU ZARANOU	
DS ABENGOUROU	DISPENSARE DIOULAKRO	
DS ABENGOUROU	MATERNITE CAFETOU	
DS ABENGOUROU	PIM ABENGOUROU	
DS ABENGOUROU	PMI ABENGOUROU	
DS ABENGOUROU	SSSU ABENGOUROU	

Figure 77. Select Sites

- 4 Click Select all to add all sites to the quantification

Or, select the individual sites to add to the quantification and click Select

Or, click Cancel button to close form and return to Working Area screen.

Next, you can indicate which categories of testing services are available at each site, and which supplies each site requires. In some countries, all services will be available at every site, and all supplies will be required. However, in other countries, certain services will only be available at some sites and not at others. For each site,

- 5 Double click on the  or  to indicate whether or not a service is available.

- 6** Repeat steps 2 through 5 for each category.
- 7** Click Save button to save working area.
- 8** Click Next button to move to next working area.

## Initial Patient Data

Before forecasting patient demand during the upcoming time period, a baseline patient level in “Time Zero” using known numbers of patients on treatment and on pre-treatment at each site must be set.



Figure 78. *Initial Patient Data*

There are two ways of setting a baseline:

- Option 1: Enter patient numbers for each site
- Option 2: Enter patient numbers for a different month that is available and apply the same percentages to sites for the current forecast

It is recommended to use Option 1 if site level data from the last month prior to the start of the forecast is available; otherwise Option 2 allows the use of the most recent data available.

- 1** Click on the checkbox to 1-Enter patient numbers for each site

Or, Click on the checkbox to 2-Enter patient numbers for a different month that is available and apply the same percentages to sites for the current forecast.

- 2** Click Save button to save working area.
- 3** Click Next button to move to next working area.

## Current Patient Numbers by Site

Depending on which option was chosen on the previous Working Area screen, the appropriate Current Patient Numbers by Site screen will display.

### OPTION 1

If Option 1 was chosen on the previous screen, this screen will appear (Figure 79). For each site, the number of patients on treatment and pre-treatment for “Time Zero” needs to be entered.

Category/Region	ART Site	Patients On Treatment	Patients On Pre-Treatment
DS ABENGOUROU	BETHESDA ABENGOUROU	0	0
DS ABENGOUROU	CAT ABENGOUROU	0	0
DS ABENGOUROU	CHR ABENGOUROU	0	0
DS ABENGOUROU	CSU AFFALIKRO	0	0
DS ABENGOUROU	CSU AMELEKIA	0	0
DS ABENGOUROU	CSU APPROMPRON	0	0
DS ABENGOUROU	CSU BETTIE	0	0
DS ABENGOUROU	CSU NIABLE	0	0
DS ABENGOUROU	CSU SANKADIOKRO	0	0
DS ABENGOUROU	CSU YAKASSE-FEYASSE	0	0
DS ABENGOUROU	CSU ZARANDU	0	0
DS ABENGOUROU	DISPENSARE DIOULAKRO	0	0
DS ABENGOUROU	MATERNITE CAFETOU	0	0
DS ABENGOUROU	PIM ABENGOUROU	0	0
DS ABENGOUROU	PMI ABENGOUROU	0	0

Figure 79. Option 1: Current Patient Numbers by Site

- 1 Double click in the Patients on Treatment cell.
- 2 Type in the number of patients on treatment
- 3 Double Click in the Patients on Pre-Treatment cell.
- 4 Type in the number of patients on pre-treatment.
- 5 Repeat steps 20.1 to 20.4 for each site on the list.

When all sites have been entered –

- 6 Click Save button to save working area
- 7 Click Next button to move to next working area.

#### OPTION 2

If Option 2 was chosen on the previous screen, this screen will appear (Figure 80). The On Treatment and On Pre-Treatment fields indicate the total number of patients for the last month prior to the start of the forecast.

Category/Region	ART Site	Patients On Treatment	Patients On Pre-Treatment	% Treatm...	% Pre-Trea...	New Targe...	New Targ...
DS ABENGOUROU	BETHESDA ABENGOUROU	0	0	0	0	0	0
DS ABENGOUROU	CAT ABENGOUROU	0	0	0	0	0	0
DS ABENGOUROU	CHR ABENGOUROU	0	0	0	0	0	0
DS ABENGOUROU	CSU AFFALIKRO	0	0	0	0	0	0
DS ABENGOUROU	CSU AMELEKIA	0	0	0	0	0	0
DS ABENGOUROU	CSU APPROMPRON	0	0	0	0	0	0
DS ABENGOUROU	CSU BETTIE	0	0	0	0	0	0
DS ABENGOUROU	CSU NIABLE	0	0	0	0	0	0
DS ABENGOUROU	CSU SANKADIOKRO	0	0	0	0	0	0

Figure 80. Option 2: Current Patient Numbers by Site

- 1 Type in the current total number of patients On Treatment.
- 2 Type in the current total number of patients On Pre-Treatment.

For each site, the number of patients on treatment and pre-treatment for the most recent month needs to be entered.

- 3 Double click in the Patients on Treatment cell.
- 4 Type in the number of patients on treatment
- 5 Double Click in the Patients on Pre-Treatment cell.
- 6 Type in the number of patients on pre-treatment.
- 7 Repeat steps 3 to 6 for each site on the list.

When all sites have been entered, the model will automatically use the percent allocation of the patients across sites from the older data to estimate the number of patients at each site for “Time Zero.”

- 8** Click Save button to save working area
- 9** Click Next button to move to next working area.

## Patient Treatment Target

In preparation for entering the patient treatment targets, users must select whether to set targets according to the number of patients that are ON TREATMENT at the end of the forecast period, or the number of patients that have EVER STARTED ON TREATMENT by the end of the forecast period. Different Countries often track the number of patients on treatment in different ways, depending on what data is available. Due to high levels of attrition and death, it is often difficult to estimate how many patients remain on treatment at any given time. The number of patients that have ever started on treatment will not take into account attrition and death, whereas the number of patients actually on treatment will. In order to reach a given target of patients on treatment, the system must add enough patients over the course of the forecast period to outweigh the patients lost due to attrition and death.



Figure 81. Patient Treatment Target

And so, there are two ways of entering a target:

- Option 1: Patients actually on treatment
  - Option 2: Patients EVER STARTED on treatment
- 1** Click on the checkbox to 1-Patients actually ON treatment

Or, Click on the checkbox to 2-Patients EVER STARTED on Treatment.

- 2** Click Save button to save working area.
- 3** Click Next button to move to next working area. If Option 1 is chosen, continue with step 1.

## Patient Ever Started on Treatment

If Option 2 is chosen, the Patient Ever Started On Treatment screen will appear in order to select a method for establishing a baseline of how many patients were EVER STARTED ON TREATMENT before the forecast time period (Figure 82).



Figure 82. Patient Ever Started on Treatment

There are three ways of setting a baseline:

- Option 1: Enter EVER STARTED patient numbers for each site
- Option 2: Enter data for a different month that is available and apply the same percentages to sites for the current forecast
- Option 3: Enter on national number for patients EVER STARTED. Apply same percentages used for current patients to each of the sites.

It is recommended to use Option 1 if site level data from the last month prior to the start of the forecast is available; Option 2 allows the use of the most recent data available if site level data from the last month prior to the start of the forecast is unavailable; otherwise Option 3 should be used if there is no data on how many patients have EVER STARTED ON TREATMENT at each site prior to the start of the forecast period.

- 1 Click on the checkbox to 1- Enter EVER STARTED patient numbers for each site

Or, Click on the checkbox to 2- Enter data for a different month that is available and apply the same percentages to sites for the current forecast.

Or, Click on the checkbox to 3- Enter on national number for patients EVER STARTED. Apply same percentages used for current patients to each of the sites.

- 2 Click Save button to save working area.
- 3 Click Next button to move to next working area.

## EVER STARTED Patient Numbers by Site

Depending on which option was chosen on the previous Working Area screen, the appropriate Current Patient Numbers by Site screen will display.

### OPTION 1

If Option 1 was chosen on the previous screen, this screen will appear (Figure 83). For each site, the number of EVER STARTED patients on treatment and pre-treatment for “Time Zero” needs to be entered.



Category/Region	ART Site	Ever Started Patients On Treatment	Ever Started Patients On Pre-Treatment
DS ABENGOUROU	BETHESDA ABENGOUROU	0	0
DS ABENGOUROU	CAT ABENGOUROU	0	0
DS ABENGOUROU	CHR ABENGOUROU	0	0
DS ABENGOUROU	CSU AFFALIKRO	0	0
DS ABENGOUROU	CSU AMELEKIA	0	0
DS ABENGOUROU	CSU APPROMPRON	0	0
DS ABENGOUROU	CSU BETTIE	0	0
DS ABENGOUROU	CSU NIABLE	0	0
DS ABENGOUROU	CSU SANKADIOKRO	0	0
DS ABENGOUROU	CSU YAKASSE-FEYASSE	0	0
DS ABENGOUROU	CSU ZARANOU	0	0
DS ABENGOUROU	DISPENSARE DIOULAKRO	0	0
DS ABENGOUROU	MATERNITE CAFETOU	0	0
DS ABENGOUROU	PIM ABENGOUROU	0	0
DS ABENGOUROU	PMI ABENGOUROU	0	0
DS ABENGOUROU	SSSU ABENGOUROU	0	0
DS ABOBO - EST	CAT ABOBO	0	0

**EVER STARTED Patient Numbers by Site**

Enter Patient Data from the Last Month Prior to the Start of the Forecast

Patients on Treatment and Pre-Treatment In the two columns of white cells, enter the number of patients on treatment and on pre-treatment at each site for 'Time Zero.'

Figure 83. Option 1: EVER STARTED Patient Numbers by Site

- 1 Double click in the Ever Started Patients on Treatment cell.
- 2 Type in the number of ever started patients on treatment
- 3 Double Click in the Ever Started Patients on Pre-Treatment cell.
- 4 Type in the number of ever started patients on pre-treatment.
- 5 Repeat steps 22.1 to 27.4 for each site on the list.

When all sites have been entered –

- 6 Click Save button to save working area
- 7 Click Next button to move to next working area.

## OPTION 2

If Option 2 was chosen on the previous screen, this screen will appear (Figure 84). The On Treatment and On Pre-Treatment fields indicate the total number of ever started patients for the last month prior to the start of the forecast.

Category/Region	ART Site	Ever Started On Tre...	EVER STARTED On Pre...	% Treat	% Pre-T	New Ta...
DS ABENGOUROU	BETHESDA ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	CAT ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	CHR ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	CSU AFFALIKRO	0	0	0	0	0
DS ABENGOUROU	CSU AMELEKIA	0	0	0	0	0
DS ABENGOUROU	CSU APPROMPRON	0	0	0	0	0
DS ABENGOUROU	CSU BETTIE	0	0	0	0	0
DS ABENGOUROU	CSU NIABLE	0	0	0	0	0
DS ABENGOUROU	CSU SANKADIOKRO	0	0	0	0	0
DS ABENGOUROU	CSU YAKASSE-FEYASSE	0	0	0	0	0
DS ABENGOUROU	CSU ZARANOU	0	0	0	0	0
DS ABENGOUROU	DISPENSARE DIOULAKRO	0	0	0	0	0
DS ABENGOUROU	MATERNITE CAFETOU	0	0	0	0	0
DS ABENGOUROU	PIM ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	PMI ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	SSSU ABENGOUROU	0	0	0	0	0
DS ABOBO - EST	CAT ABOBO	0	0	0	0	0
DS ABOBO - EST	CES ABOBOTE	0	0	0	0	0
DS ABOBO - EST	CSS ANGE GARDIEN	0	0	0	0	0
DS ABOBO - EST	CSU BANCO SUD	0	0	0	0	0
DS ABOBO - EST	CSU COM ABOBOTE	0	0	0	0	0
DS ABOBO - EST	CSU COM KENNEDY CLOUET	0	0	0	0	0
DS ABOBO - EST	FSU COM ABOBO AVOCATIER	0	0	0	0	0
DS ABOBO - EST	FSU COM ABOBO BAOULE	0	0	0	0	0
DS ABOBO - EST	FSU COM ABOBO BC	0	0	0	0	0

**EVER STARTED Patient Numbers by Site**

**Total Patients**

You must enter the total number of patients on treatment and in pre-treatment across all sites for the last month prior to the start of the forecast

**Patients on Treatment and Pre-Treatment**

In the two columns, enter the number of patients on treatment and on pre-treatment at each site for the most recent month you have available. The model will then automatically use the percent allocation of patients across sites from the older data to estimate the number of patients at each site for 'Time Zero.'

Figure 84. Option 2: EVER STARTED Patient Numbers by Site

- 1 Type in the current total number of ever started patients On Treatment.
- 2 Type in the current total number of ever started patients On Pre-Treatment.

For each site, the number of ever started patients on treatment and pre-treatment for the most recent month needs to be entered.

- 3 Double click in the Ever Started on Treatment cell.
- 4 Type in the number of ever started patients on treatment
- 5 Double Click in the Ever Started on Pre-Treatment cell.
- 6 Type in the number of ever started patients on pre-treatment.
- 7 Repeat steps 3 to 6 for each site on the list.

When all sites have been entered, the model will automatically use the percent allocation of the patients across sites from the older data to estimate the number of patients at each site for “Time Zero.”

- 8 Click Save button to save working area
- 9 Click Next button to move to next working area.

### OPTION 3

If Option 3 was chosen on the previous screen, this screen will appear (Figure 85). The On Treatment and On Pre-Treatment fields indicate the total number of ever started patients for the last month prior to the start of the forecast.

The screenshot shows the 'Morbidity Quantification' tool interface. The left sidebar contains the following sections:

- Forecast Description:**
  - Title of Quantification: Morbidity Quantification
  - Date of quantification: Friday, December 27, 2013
  - Description, purpose and key assumption: This quantification is for demonstration purposes only.
- Forecast parameters:**
  - Forecast start date: 1/ 1/2014
  - Budget Start: January
  - Budget End: December
  - Type of Algorithm: Serial (selected), Parallel

The main area of the tool features a table with the following columns: Category/Region, A/T Site, % Treatment, % Pre-Treatment, New Target Treat, and New Target Pre-Treatment. Above the table, there are fields for 'Total number of patients' and 'How many patients EVER STARTED treatment by', with sub-fields for 'On Treatment' and 'On Pre-Treatment'. The table is currently empty, and the 'On Treatment' and 'On Pre-Treatment' fields are also empty.

Figure 85. Option 3: EVER STARTED Patient Numbers by Site

- 1 Type in the current total number of ever started patients On Treatment.
- 2 Type in the current total number of ever started patients On Pre-Treatment.

The model will automatically use the percent allocation of the patients across sites from the current data to estimate the number of patients at each site for “Time Zero.”



- 3** Click Save button to save working area
- 4** Click Next button to move to next working area.

## ART Patient Target

Next, the preferred method of setting ART patient targets for the selected site groups in the quantification is selected.



Figure 86. ART Patient Target

There are four options depending on how much data is available:

- Option 1: Enter National Targets – Simply apply equal linear growth to all sites
- Option 2: Enter National Targets – But grow some sites at a different rate
- Option 3: Enter National Targets – But enter specific targets for some sites
- Option 4: Enter specific targets for each site

It is recommended to select Option 1 if you can only enter a national target and you plan to assume that all sites will grow at the same rate to hit this target; Option 2 if you want to enter one national target but want the flexibility to manually set different growth rates for different sites; Option 3 if you want to enter one national target but want the flexibility to manually set different absolute targets for different sites; otherwise, Option 4 if you want to manually set a different target for every site without entering a national target.

- 1** Click on the checkbox to 1 – Enter National Targets – Simply apply equal linear growth to all sites

Or, Click on the checkbox to 2 – Enter National Targets – But grow some sites at a different rate

Or, Click on the checkbox to 3 – Enter National Targets – But enter specific targets for some sites

Or, Click on the checkbox to 4 – Enter specific targets for each site

- 2** Click Save button to save working area.
- 3** Click Next button to move to next working area. If Option 4 is chosen, continue with step 1.

## Site Target Calculator

Depending on which option was chosen on the previous Working Area screen, the appropriate Site Target Calculator screen will display.

### OPTION 1

If Option 1 was chosen on the previous screen, this screen will appear (Figure 87). For each month, the National Treatment Target needs to be entered.

The screenshot shows the 'Site Target Calculator' window. At the top, there are buttons for 'Home', 'Back', 'Site Target Calculator', and 'Next'. Below this is a section titled 'National Treatment Target' containing a table with 13 columns: December, January, February, March, April, May, June, July, August, September, October, November, and December, followed by an 'Annual Growth' column. The table is populated with the value '1367081' for all months and '0%' for Annual Growth. Below the table, there is a field labeled '% of Patients Who Are Children' with the value '0'. To the right of this field are two buttons: 'Apply Pediatric % to All Sites' and 'Calculate linear scale-up between December 2014 December 2015'.

December	January	February	March	April	May	June	July	August	September	October	November	December	Annual Growth
1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	0%

% of Patients Who Are Children: 0

Buttons: Apply Pediatric % to All Sites, Calculate linear scale-up between December 2014 December 2015

Figure 87. Option 1: Site Target Calculator

1. Enter the number of patients expected to be on treatment by the end of the forecast period in the white cell
2. Populate the optional green cells by manually inputting monthly national targets

Or, Click on the Calculate Linear Scale Up button to assume linear growth.

To allow for demographic differences between adults and children, the percentage of patients on treatment who are children for the group of sites in the quantification should be entered.

3. Enter the % of Patients who are children.
4. Click Apply Pediatric % to All Sites.
5. Click Save button to save working area
6. Click Next button to move to next working area.

### OPTION 2

If Option 2 was chosen on the previous screen, this screen will appear (Figure 88). For each month, the National Treatment Target needs to be entered.

Home Back Site Target Calculator Next

National Treatment Target

December	January	February	March	April	May	June	July	August	September	October	November	December	Annual Growth
1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	1367081	0%

% of Patients Who Are Children: 0

Apply Pediatric % to All Sites

Calculate linear scale-up between December 2014 December 2015

Category/Region	ART Site	% pediatrics	Dec...	Jan...	Feb...	Mar...	April	May	June	July	August	Sept...	Oct...	Nov...	Dec...	Growth Target	Apply Growth
DS ABENGOURDI	BETHESDA AB 0		2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CAT ABENGOL 0		63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CHR ABENGOL 0		2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CSU AFFALIKR 0		32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CSU AMELEKI 0		64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CSU APPROMF 0		86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CSU BETTIE 0		56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CSU NIABLE 0		45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	0	<input checked="" type="checkbox"/>
DS ABENGOURDI	CSU SANKADI 0		4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	0	<input checked="" type="checkbox"/>

Figure 88. Option 2: Site Target Calculator

- 1 Enter the number of patients expected to be on treatment by the end of the forecast period in the white cell
- 2 Populate the optional green cells by manually inputting monthly national targets

Or, Click on the Calculate Linear Scale Up button to assume linear growth.

To allow for demographic differences between adults and children, the percentage of patients on treatment who are children for the group of sites in the quantification should be entered.

- 3 Enter the % of Patients who are children.
- 4 Click Apply Pediatric % to All Sites.

The lower part of the screen lists the ART sites in the quantification and displays the targets for each month in the forecast including the % of patients who are children at the site and the growth rate at the site.

The % pediatrics can be entered separately if the percentage is different than the national average for that specific site.

- 5 Enter % Pediatrics for each site where pediatric data at site level is available. (For a site where additional data is unavailable, assume the same percentage of pediatrics as the overall population.)

The growth target rate can be entered separately if the rate is different than the national average for that specific site.

- 6 Enter Growth Target for each site where the % growth rate at site level is available. (For a site where additional data is unavailable, assume the same growth rate as the national average.)
- 7 Check the box next to Growth Target field to apply the site specific target.
- 8 Click Save button to save working area
- 9 Click Next button to move to next working area.

**OPTION 3**

If Option 3 was chosen on the previous screen, this screen will appear (Figure 89). For each month, the National Treatment Target needs to be entered.

The screenshot shows the 'Site Target Calculator' window. At the top, there are buttons for 'Home', 'Back', 'Site Target Calculator', and 'Next'. Below these is the 'National Treatment Target' section, which includes a grid for entering monthly targets (December to December) and an 'Annual Growth' field set to 0%. There is also a '% of Patients Who Are Children' field set to 0 and a button 'Apply Pediatric % to All Sites'. A button 'Calculate linear scale-up between December 2014 December 2015' is also present. Below this is a table listing ART sites with their targets for each month and a 'Growth Target' field.

Category/Region	ART Site	% pediatrics	Dec...	Jan...	Feb...	Mar...	April	May	June	July	August	Sept...	Oct...	Nov...	Dec...	Growth Target	Apply Growth
DS ABENGOUROI	BETHESDA AB 0		2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	2346	0	<input type="checkbox"/>
DS ABENGOUROI	CAT ABENGOL 0		63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	0	<input type="checkbox"/>
DS ABENGOUROI	CHR ABENGOL 0		2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	0	<input type="checkbox"/>
DS ABENGOUROI	CSU AFFALIKR 0		32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	0	<input type="checkbox"/>
DS ABENGOUROI	CSU AMELEKI 0		64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	0	<input type="checkbox"/>
DS ABENGOUROI	CSU APPROMF 0		86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	86234	0	<input type="checkbox"/>
DS ABENGOUROI	CSU BETTIE 0		56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	56235	0	<input type="checkbox"/>
DS ABENGOUROI	CSU NIABLE 0		45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	45234	0	<input type="checkbox"/>
DS ABENGOUROI	CSU SANKADI 0		4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	4534	0	<input type="checkbox"/>

Figure 89. Option 3: Site Target Calculator

- 1 Enter the number of patients expected to be on treatment by the end of the forecast period in the white cell
- 2 Populate the optional green cells by manually inputting monthly national targets

Or, Click on the Calculate Linear Scale Up button to assume linear growth.

To allow for demographic differences between adults and children, the percentage of patients on treatment who are children for the group of sites in the quantification should be entered.

- 3 Enter the % of Patients who are children.
- 4 Click Apply Pediatric % to All Sites.

The lower part of the screen lists the ART sites in the quantification and displays the targets for each month in the forecast including the % of patients who are children at the site and the growth rate at the site.

The % pediatrics can be entered separately if the percentage is different than the national average for that specific site.


- 5 Enter % Pediatrics for each site where pediatric data at site level is available. (For a site where additional data is unavailable, assume the same percentage of pediatrics as the overall population.)

The actual target level can be entered separately if the target is different than the national average for that specific site.

- 6 Enter Growth Target for each site where the actual target level at site level is available. (For a site where additional data is unavailable, assume the same target as the calculated value.)
- 7 Check the box next to Growth Target field to apply the site specific target.
- 8 Click Save button to save working area
- 9 Click Next button to move to next working area.

## Site Specific Treatment Patient Numbers

The next screen (Figure 90) allows the pre-populated targets by month and by site from the Site Target Calculator screen to be overridden. If Option 4 was chosen on the ART Patient Target screen, then this is the screen where the detailed targets by month and by site can be entered.



Category/Region	ART Site	% pedia...	Decem...	January	February	March	April	May	June	July	August
DS ABENGOUROU	BETHESDA ABENGOUROU	0	2346	2346	2346	2346	2346	2346	2346	2346	2346
DS ABENGOUROU	CAT ABENGOUROU	0	63347	63347	63347	63347	63347	63347	63347	63347	63347
DS ABENGOUROU	CHR ABENGOUROU	0	2344	2344	2344	2344	2344	2344	2344	2344	2344
DS ABENGOUROU	CSU AFFALIKRO	0	32445	32445	32445	32445	32445	32445	32445	32445	32445
DS ABENGOUROU	CSU AMELEKIA	0	64324	64324	64324	64324	64324	64324	64324	64324	64324
DS ABENGOUROU	CSU APPROMPRON	0	86234	86234	86234	86234	86234	86234	86234	86234	86234
DS ABENGOUROU	CSU BETTIE	0	56235	56235	56235	56235	56235	56235	56235	56235	56235
DS ABENGOUROU	CSU NIABLE	0	45234	45234	45234	45234	45234	45234	45234	45234	45234
DS ABENGOUROU	CSU SANKADIOKRO	0	4534	4534	4534	4534	4534	4534	4534	4534	4534
DS ABENGOUROU	CSU YAKASSE-FEYASE	0	67563	67563	67563	67563	67563	67563	67563	67563	67563
DS ABENGOUROU	CSU ZARANDOU	0	7566	7566	7566	7566	7566	7566	7566	7566	7566
DS ABENGOUROU	DISPENSARE DIOULA	0	4556	4556	4556	4556	4556	4556	4556	4556	4556
DS ABENGOUROU	MATERNITE CAFETOU	0	654	654	654	654	654	654	654	654	654
DS ABENGOUROU	PIM ABENGOUROU	0	765	765	765	765	765	765	765	765	765
DS ABENGOUROU	PMI ABENGOUROU	0	4565	4565	4565	4565	4565	4565	4565	4565	4565
DS ABENGOUROU	SSSU ABENGOUROU	0	676	676	676	676	676	676	676	676	676
DS ABOBO - EST	CAT ABOBO	0	456	456	456	456	456	456	456	456	456
DS ABOBO - EST	CES ABOBO	0	86	86	86	86	86	86	86	86	86
DS ABOBO - EST	CSS ANGE GARDIEN	0	65464	65464	65464	65464	65464	65464	65464	65464	65464

Figure 90. Site Specific Treatment Patient Numbers

- 1 Double click in the % pediatrics column.
- 2 Type in the % of patients who are children at the site.
- 3 Repeat steps 32 and 33 for each site in the quantification.
- 4 Double click in the target patient number field for the appropriate month.
- 5 Type in the target patient number for that month at that site.
- 6 Repeat for each month in the forecast and each site in the quantification.
- 7 Click Save button to save working area.
- 8 Click Next button to move to next working area.

## Patient Pre-Treatment Targets

Next, the preferred method of setting ART patient targets for the selected site groups in the quantification is selected.

Figure 91. Patient Pre-Treatment Targets

There are four options depending on how much data is available:

- Option 1: Enter National Pre-Treatment Targets
- Option 2: Enter National Pre-Treatment Targets – But grow some sites at a different rate
- Option 3: Enter National Pre-Treatment Targets – But enter specific targets for some sites
- Option 4: Enter specific Pre-Treatment targets for each site
- Option 5: No Pre-Treatment targets are available. Use Testing Efficiency to estimate.

It is recommended to select Option 1 if you can only enter a national target and you plan to assume that all sites will grow at the same rate to hit this target; Option 2 if you want to enter one national target but want the flexibility to manually set different growth rates for different sites; Option 3 if you want to enter one national target but want the flexibility to manually set different absolute targets for different sites; Option 4 if you want to manually set a different target for every site without entering a national target; otherwise Option 5 if you do not want to input any pre-treatment targets and will use Testing Efficiency instead.

- 1 Click on the checkbox to 1 – Enter National Pre-Treatment Targets – Simply apply equal linear growth to all sites

Or, Click on the checkbox to 2 – Enter National Pre-Treatment Targets – But grow some sites at a different rate

Or, Click on the checkbox to 3 – Enter National Pre-Treatment Targets – But enter specific targets for some sites

Or, Click on the checkbox to 4 – Enter specific Pre-Treatment targets for each site

Or, Click on the checkbox to 5 – No Pre-Treatment Targets are available. Use Testing Efficiency to estimate.

- 2 Click Save button to save working area.
- 3 Click Next button to move to next working area. . If Option 4 is chosen, continue with step 1. If Option 5 is chosen, continue with step 1.

## Site Pre-Treatment Target Calculator

Depending on which option was chosen on the previous Working Area screen, the appropriate Site Pre-Treatment Target Calculator screen.

### OPTION 1

If Option 1 was chosen on the previous screen, this screen will appear (Figure 92). For each month, the National Pre-Treatment Target needs to be entered.

Figure 92. Option 1: Site Pre-Treatment Target Calculator

- 1 Enter the number of patients expected to be on pre-treatment by the end of the forecast period in the white cell
- 2 Populate the optional green cells by manually inputting monthly national targets

Or, Click on the Calculate Linear Scale Up button to assume linear growth.

To allow for demographic differences between adults and children, the percentage of patients on pre-treatment who are children for the group of sites in the quantification should be entered.

- 3 Enter the % of Patients who are children.
- 4 Click Apply Pediatric % to All Sites.
- 5 Click Save button to save working area
- 6 Click Next button to move to next working area.

### OPTION 2

If Option 2 was chosen on the previous screen, this screen will appear (Figure 93). For each month, the National Pre-Treatment Target needs to be entered.

Figure 93. Option 2: Site Pre-Treatment Target Calculator

- 1 Enter the number of patients expected to be on pre-treatment by the end of the forecast period in the white cell
- 2 Populate the optional green cells by manually inputting monthly national targets



Or, Click on the Calculate Linear Scale Up button to assume linear growth.

To allow for demographic differences between adults and children, the percentage of patients on pre-treatment who are children for the group of sites in the quantification should be entered.

**3** Enter the % of Patients who are children.

**4** Click Apply Pediatric % to All Sites.

The lower part of the screen lists the ART sites in the quantification and displays the targets for each month in the forecast including the % of patients who are children at the site and the growth rate at the site.

The % pediatrics can be entered separately if the percentage is different than the national average for that specific site.

**5** Enter % Pediatrics for each site where pediatric data at site level is available. (For a site where additional data is unavailable, assume the same percentage of pediatrics as the overall population.)

The growth target rate can be entered separately if the rate is different than the national average for that specific site.

**6** Enter Growth Target for each site where the % growth rate at site level is available. (For a site where additional data is unavailable, assume the same growth rate as the national average.)

**7** Check the box next to Growth Target field to apply the site specific target.

**8** Click Save button to save working area

**9** Click Next button to move to next working area.

### OPTION 3

If Option 3 was chosen on the previous screen, this screen will appear (Figure 94). For each month, the National Pre-Treatment Target needs to be entered.

The screenshot shows the 'Site Target Calculator' window. At the top, there are buttons for 'Home', 'Back', 'Site Target Calculator', and 'Next'. Below these is the 'National Treatment Target' section, which includes a table of monthly targets (December to December) and an 'Annual Growth' field set to 0%. Below this is a section for '% of Patients Who Are Children' with a value of 0 and a button 'Apply Pediatric % to All Sites'. To the right of this is a button 'Calculate linear scale-up between December 2014 December 2015'. Below these sections is a table with columns for 'Category/Region', 'ART Site', '% pediatrics', and months from Dec to Dec, followed by 'Growth Target' and 'Apply Growth' checkboxes. The table lists four ART sites: D5 ABENGOURDI BETHESDA AB 0, D5 ABENGOURDI CAT ABENGOL 0, D5 ABENGOURDI CHR ABENGOL 0, and D5 ABENGOURDI CSU AFFALKR 0. Each site has a row of monthly targets and a growth target of 0.

Category/Region	ART Site	% pediatrics	Dec	Jan	Feb	Mar	April	May	June	July	August	Sept	Oct	Nov	Dec	Growth Target	Apply Growth
D5 ABENGOURDI BETHESDA AB 0			2345	2345	2345	2345	2345	2345	2345	2345	2345	2345	2345	2345	2345	0	<input type="checkbox"/>
D5 ABENGOURDI CAT ABENGOL 0			63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	63347	0	<input type="checkbox"/>
D5 ABENGOURDI CHR ABENGOL 0			2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	2344	0	<input type="checkbox"/>
D5 ABENGOURDI CSU AFFALKR 0			32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	32445	0	<input type="checkbox"/>
D5 ABENGOURDI CSU AMELEKSI 0			64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	64324	0	<input type="checkbox"/>

Figure 94. Option 3: Site Pre-Treatment Target Calculator

**1** Enter the number of patients expected to be on pre-treatment by the end of the forecast period in the white cell

**2** Populate the optional green cells by manually inputting monthly national targets

Or, Click on the Calculate Linear Scale Up button to assume linear growth.

To allow for demographic differences between adults and children, the percentage of patients on pre-treatment who are children for the group of sites in the quantification should be entered.



**3** Enter the % of Patients who are children.

**4** Click Apply Pediatric % to All Sites.

The lower part of the screen lists the ART sites in the quantification and displays the targets for each month in the forecast including the % of patients who are children at the site and the growth rate at the site.

The % pediatrics can be entered separately if the percentage is different than the national average for that specific site.

**5** Enter % Pediatrics for each site where pediatric data at site level is available. (For a site where additional data is unavailable, assume the same percentage of pediatrics as the overall population.)

The actual target level can be entered separately if the target is different than the national average for that specific site.

**6** Enter Growth Target for each site where the actual target level at site level is available. (For a site where additional data is unavailable, assume the same target as the calculated value.)

**7** Check the box next to Growth Target field to apply the site specific target.

**8** Click Save button to save working area

**9** Click Next button to move to next working area.

## Site Specific Pre-Treatment Patient Numbers

The next screen (Figure 95) allows the pre-populated targets by month and by site from the Site Pre-Treatment Target Calculator screen to be overridden. If Option 4 was chosen on the Patient Pre-Treatment Target screen, then this is the screen where the detailed targets by month and by site can be entered.

Category/Region	ART Site	% pediatric	Decem	January	February	March	April	May	June	July	August
DS ABENGOUROU	BETHESDA ABENGOUROU 0		656	656	656	656	656	656	656	656	656
DS ABENGOUROU	CAT ABENGOUROU 0		3454	3454	3454	3454	3454	3454	3454	3454	3454
DS ABENGOUROU	CHI ABENGOUROU 0		7634	7634	7634	7634	7634	7634	7634	7634	7634
DS ABENGOUROU	CSU AFFALIKPO 0		863	863	863	863	863	863	863	863	863
DS ABENGOUROU	CSU AMELEKIA 0		44563	44563	44563	44563	44563	44563	44563	44563	44563
DS ABENGOUROU	CSU APPROMPRON 0		34563	34563	34563	34563	34563	34563	34563	34563	34563
DS ABENGOUROU	CSU BETTIE 0		6546	6546	6546	6546	6546	6546	6546	6546	6546
DS ABENGOUROU	CSU NAMIE 0		7356	7356	7356	7356	7356	7356	7356	7356	7356
DS ABENGOUROU	CSU SANKADIKRO 0		834	834	834	834	834	834	834	834	834
DS ABENGOUROU	CSU YAKASSE FEYASE 0		8856	8856	8856	8856	8856	8856	8856	8856	8856
DS ABENGOUROU	CSU ZARANDU 0		67	67	67	67	67	67	67	67	67
DS ABENGOUROU	DISPENSARIE DIOULA 0		53456	53456	53456	53456	53456	53456	53456	53456	53456
DS ABENGOUROU	MATERINITE CAFETOU 0		43345	43345	43345	43345	43345	43345	43345	43345	43345
DS ABENGOUROU	PIM ABENGOUROU 0		4634	4634	4634	4634	4634	4634	4634	4634	4634
DS ABENGOUROU	PMI ABENGOUROU 0		65634	65634	65634	65634	65634	65634	65634	65634	65634

Figure 95. Site Specific Pre-Treatment Patient Numbers

**1** Double click in the % pediatrics column.

**2** Type in the % of patients who are children at the site.

**3** Repeat steps 32 and 33 for each site in the quantification.

**4** Double click in the target pre-treatment patient number field for the appropriate month.

**5** Type in the target pre-treatment patient number for that month at that site.

**6** Repeat for each month in the forecast and each site in the quantification.

**7** Click Save button to save working area.

**8** Click Next button to move to next working area. Continue with step 1

## Testing Efficiency

If you selected Option 1 through 4 on the previous screen, this screen will be skipped.

Testing efficiency is the percentage of HIV-positive patients who qualify for initiating ART. It can be unique to a specific country or site, and can be used to estimate how many patients will be receiving pre-treatment care assuming a given level of patients receiving treatment. This estimate will also take into consideration overall prevalence rates of HIV in the country or sites in question. Although this method can produce a useful estimate, it is only recommended if you do not have more substantial data at the site level.

Category/Region	ART Sites	A.T.E.	P.T.E.	% Pediatric
D5 ABENGOUROU	BETHESDA ABENGOUROU	0	0	0
D5 ABENGOUROU	CAT ABENGOUROU	0	0	0
D5 ABENGOUROU	CHR ABENGOUROU	0	0	0
D5 ABENGOUROU	CSU AFFALIKRO	0	0	0
D5 ABENGOUROU	CSU AMELEKA	0	0	0
D5 ABENGOUROU	CSU APPROMPRON	0	0	0
D5 ABENGOUROU	CSU BETTIE	0	0	0

Figure 96. Testing Efficiency

- 1 Type in the default Adult Testing Efficiency (A.T.E.) value.
- 2 Type in the default Pediatric Testing Efficiency (P.T.E.) value.
- 3 Type in the default % Pediatric (Pre-Existing Patients).

The default values can be loaded as the default value for all sites. Each site can then be edited as needed.

- 4 Click the Apply Default Values to All Sites button.

For each site in the quantification –

- 5 Update the A.T.E., P.T.E., or % Pediatric fields as needed.
- 6 Click Save button to save working area.
- 7 Click Next button to move to next working area.

## Testing Information

Category/Region	ART Sites	A.T.P.	P.T.P.	A.D.D.W.F.	P.D.D.W.F.	D.T.F.W.R.
DS ABENGOUROU	BETHESDA ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	CAT ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	CHR ABENGOUROU	0	0	0	0	0
DS ABENGOUROU	CSU AFFALIKRO	0	0	0	0	0
DS ABENGOUROU	CSU AMELEKIA	0	0	0	0	0
DS ABENGOUROU	CSU APPROMPRON	0	0	0	0	0
DS ABENGOUROU	CSU BETTIE	0	0	0	0	0

Figure 97. Testing Information

The % Adult Testing Population HIV Positive field (A.T.P.) is the percentage of adults test for HIV that receives a positive test result.

- 1 Type in the default % Adult Testing Population HIV Positive value.

The % Pediatric Testing Population HIV Positive field (P.T.P.) is the percentage of children tested for HIV that receive a positive test result.

- 2 Type in the default % Pediatric Testing Population HIV Positive value.

The % of HIV+ Adult diagnoses to depart w/out follow-up field (A.D.D.W.F.) is the percentage of adults with positive HIV test results that fail to follow up for care after diagnosis.

- 3 Type in the default % of HIV+ Adult diagnoses to depart w/out follow-up

The % of HIV+ Pediatric diagnoses to depart w/out follow-up field (P.D.D.W.F.) is the percentage of children with positive HIV test results that fail to follow up for care after diagnosis.

- 4 Type in the default % of HIV+ Pediatric diagnoses to depart w/out follow-up

The % of HIV+ Diagnoses That Follow up Which Receive CD4 field (D.T.F.W.R.) is the percentage of HIV positive follow-up patients who are administered a CD4 test to screen for eligibility for treatment initiation.

Type in the default % of HIV+ Diagnoses That Follow up Which Receive CD4

The default values can be loaded as the default value for all sites. Each site can then be edited as needed.

- 5 Click the Apply default Values to All Sites button.

For each site in the quantification –

- 6 Update the A.T.P., P.T.P., A.D.D.W.F., P.D.D.W.F., and D.T.F.W.R. fields as needed.
- 7 Click Save button to save working area.
- 8 Click Next button to move to next working area.

## Adult Patient Behavior

Category/Region	ART Sites	AT-A.P.A.	AT-A.B.D/E.	AT-A.B.D/N.P.	APA-A.P.A.	APA-A.M.T.	APA-A.B.D/E.	APA-A.B.D/N.P.
DS ABENGOUROU	BETHESDA ABENGOUR	0	0	0	0	0	0	0
DS ABENGOUROU	CAT ABENGOUROU	0	0	0	0	0	0	0
DS ABENGOUROU	CHR ABENGOUROU	0	0	0	0	0	0	0
DS ABENGOUROU	CSU AFFALOU	0	0	0	0	0	0	0
DS ABENGOUROU	CSU AMELEKA	0	0	0	0	0	0	0
DS ABENGOUROU	CSU APPROMPRON	0	0	0	0	0	0	0
DS ABENGOUROU	CSU BETTIE	0	0	0	0	0	0	0

Figure 98. Adult Patient Behavior

For Adults in Treatment, the Annual % Patient Attrition field (AT-A.P.A.) is the percentage of adult HIV patients on treatment that terminate their treatment for any reason (including death) over the course of a year.

- 1 Type in the default Annual % Patient Attrition for Adults in treatment.

For Adults in Treatment, the Avg # Blood Draws/Existing Patient/Yr field (AT-A.B.D/E.P/Yr) is the number of blood draws that each existing adult patient on treatment receives per year. This information is critical in order to efficiently combine all applicable tests which occur at the same appointment into the same blood draw.

- 2 Type in the default Avg # Blood Draws/Existing Patient/Yr for Adults in treatment.

For Adults in Treatment, the Avg # Blood Draws/New Patient/Yr field (AT-A.B.D/N.P/Yr) is the number of blood draws that each new adult patient on treatment received per year.

- 3 Type in the default Avg # Blood Draws/New Patient/Yr for Adults in treatment.

For Adults in Pre-ART, the Annual % Patient Attrition field (APA-A.P.A.) is the percentage of adult HIV patients on pre-treatment care that terminate their care for any reason (including death).

- 4 Type in the default Annual % Patient Attrition for Adults in Pre-ART.

For Adults in Pre-ART, the Annual % Migration into Treatment field (APA-A.M.T.) is the percentage of adult HIV patients on pre-ART that migrate to ART over the course of a year.

- 5 Type in the default Annual % Migration into Treatment for Adults in Pre-ART.

For Adults in Pre-ART, the Avg # Blood Draws/Existing Patient/Yr field (APA-A.B.D/E.P/Yr) is the number of blood draws that each existing adult patient on pre-ART receives per year.

- 6 Type in the default Avg # Blood Draws/Existing Patient/Yr for Adults in Pre-ART.

For Adults in Pre-ART, the Avg # Blood Draws/New Patient/Yr field (APA -A.B.D/N.P/Yr) is the number of blood draws that each new adult patient on pre-ART received per year.

- 7** Type in the default Avg # Blood Draws/New Patient/Yr for Adults in Pre-ART.

The default values can be loaded as the default value for all sites. Each site can then be edited as needed.

- 8** Click the Apply default Values to All Sites button.

For each site in the quantification –

- 9** Update the AT-A.P.A., AT-A.B.D/E.P/Yr, AT-A.B.D/N.P/Yr, APA-A.P.A., APA-A.M.T., APA-A.B.D/E.P/Yr, and APA -A.B.D/N.P/Yr fields as needed.
- 10** 10. Click Save button to save working area.
- 11** 11. Click Next button to move to next working area.

## Pediatric Patient Behavior

Category/Region	ART Sites	PT-A.P.A.	PT-A.B.D/E.P.	PT-A.B.D/N.P.	PPA-A.P.A.	PPA-A.M.T.	PPA-A.B.D/E.	PPA-A.B.D/N.
DS ABENGOUROU	BETHESDA ABENGOUROU	0	0	0	0	0	0	0
DS ABENGOUROU	CAT ABENGOUROU	0	0	0	0	0	0	0
DS ABENGOUROU	CHR ABENGOUROU	0	0	0	0	0	0	0
DS ABENGOUROU	CSU AFFALIKRO	0	0	0	0	0	0	0
DS ABENGOUROU	CSU AMELEKIA	0	0	0	0	0	0	0
DS ABENGOUROU	CSU APPROMPRON	0	0	0	0	0	0	0
DS ABENGOUROU	CSU BETTIE	0	0	0	0	0	0	0

Figure 99. *Pediatric Patient Behavior*

For Pediatrics in Treatment, the Annual % Patient Attrition field (PT-A.P.A.) is the percentage of pediatric HIV patients on treatment that terminate their treatment for any reason (including death) over the course of a year.

- 1** Type in the default Annual % Patient Attrition for Pediatrics in Treatment.

For Pediatrics in Treatment, the Avg # Blood Draws/Existing Patient/Yr field (PT-A.B.D/E.P/Yr) is the number of blood draws that each existing pediatric patient on treatment receives per year. This information is critical in order to efficiently combine all applicable tests which occur at the same appointment into the same blood draw.

- 2** Type in the default Avg # Blood Draws/Existing Patient/Yr for Pediatrics in Treatment.

For Pediatrics in Treatment, the Avg # Blood Draws/New Patient/Yr field (PT-A.B.D/N.P/Yr) is the number of blood draws that each new pediatric patient on treatment received per year.

- 3** Type in the default Avg # Blood Draws/New Patient/Yr for Pediatrics in Treatment

For Pediatrics in Pre-ART, the Annual % Patient Attrition field (PPA-A.P.A.) is the percentage of pediatric HIV patients on pre-treatment care that terminate their care for any reason (including death).

- 4 Type in the default Annual % Patient Attrition for Pediatrics in Pre-ART.

For Pediatrics in Pre-ART, the Annual % Migration into Treatment field (PPA-A.M.T.) is the percentage of pediatric HIV patients on pre-ART that migrate to ART over the course of a year.

- 5 Type in the default Annual % Migration into Treatment for Pediatrics in Pre-ART.

For Pediatrics in Pre-ART, the Avg # Blood Draws/Existing Patient/Yr field (PPA-A.B.D/E.P/Yr) is the number of blood draws that each existing pediatric patient on pre-ART receives per year.

- 6 Type in the default Avg # Blood Draws/Existing Patient/Yr for Pediatrics in Pre-ART.

For Pediatrics in Pre-ART, the Avg # Blood Draws/New Patient/Yr field (PPA-A.B.D/N.P/Yr) is the number of blood draws that each new pediatric patient on pre-ART received per year.

- 7 Type in the default Avg # Blood Draws/New Patient/Yr for Pediatrics in Pre-ART.

The default values can be loaded as the default value for all sites. Each site can then be edited as needed.

- 8 Click the Apply default Values to All Sites button.

For each site in the quantification –

- 9 Update the PT-A.P.A., PT -A.B.D/E.P/Yr, PT -A.B.D/N.P/Yr, PPA-A.P.A., PPA-A.M.T., PPA-A.B.D/E.P/Yr, and PPA -A.B.D/N.P/Yr fields as needed.

- 10 Click Save button to save working area.

- 11 Click Next button to move to next working area.

## Rapid Test Algorithm by Site

If “Serial” algorithm was selected when defining the morbidity quantification Figure will display. For each Test in the Rapid Test Algorithm,



Figure 100. Serial Rapid Test Algorithm by Site

- 1 Type in the default % of screening, confirmatory, and tie breaker test.

The default values can be loaded as the default value for all sites. Each site can then be edited as needed.

- 2 Click the Apply default Values to All Sites button.

For each site in the quantification –

- 3 Update the % of screening, confirmatory, and tie breaker test fields as needed.
- 4 Click Save button to save working area.
- 5 Click Next button to move to next working area.

If “Parallel” algorithm was selected when defining the morbidity quantification Figure will display.

Figure 101. Parallel Rapid Test Algorithm by Site

- 6 Double click on the % value cells to enter the amount to screening and tie-breaker tests.
- 7 Repeat for each site.

## Ordering and Inventory Assumptions

Next the Ordering and Inventory Assumption screen (Figure 102) is displayed.

Figure 102. Ordering and Inventory Assumptions

The Months of Security Stock field indicates the number of months of security stock intended to include in the quantification. If you have run a national quantification in past years and have used it to effectively order sufficient quantities of supplies, there may be an existing security stock already on hand. In this case, it would not be necessary to include security stock here.

- 1 Type in the Months of Security Stock.
- 2 Repeat for each site.

The Reagents lost through wastage/leakage, etc. field indicates what percentage of reagents and other supplies are typically lost due to wastage, leakage, theft, expiration, etc. This is critical to ensure that sufficient quantities are ordered to account for this loss.

- 3 Type in the Reagents lost through wastage/leakage, etc. for each category of testing supplies.
- 4 Click Next button to move to next working area.



## Reviewing Data

The final screen of the working area is Review Your Data screen (Figure 103).

Figure 103. Review Your Data

The final step for using the morbidity methodology is to review the data that has been entered.

### PATIENT BEHAVIOR

Patient Behavior tab helps to monitor if attrition and migration percentages are above or below some specific values.

With the Review you Data screen displayed—

- 1 Click on the Patient Behavior tab. (Figure 104)

Figure 104. Patient Behavior tab

Using the Cell with error color code legend, you can determine which sites has migration rate less than a minimum and which has an attrition rate greater than a maximum value. To change the minimum and maximum value—

- 2 Double click in the % Min Migration Rate < field.
- 3 Enter the minimum migration rate.
- 4 Double click in the % Max Attrition rate > field.

- 5 Enter the maximum attrition rate.
- 6 Click Recheck button to apply the legend to the data below.

### SERVICE AND PLATFORM

Service and Platform tab helps to identify sites who doesn't have instrument for the selected testing areas or if a referral site does not have instrument for the referral service used. From this tab you can easily modify the site information.

- 1 Click Service and Platform tab.

The Service and Platform tab is displayed. (Figure 105)

Test Category	Test Select...	Test Refere...	Has Instrum...	# Instrument	Remarks
DS ABENGOUROU-> BETHESDA ABENGOUROU					
CD4	✓	✗	✗	0	Site has not instrument listed under this test category.
Chemistry	✓	✗	✗	0	Site has not instrument listed under this test category.
Hematology	✓	✗	✗	0	Site has not instrument listed under this test category.
VirallLoad	✓	✗	✗	0	Site has not instrument listed under this test category.
OtherTest	✓	✗	✗	0	Site has not instrument listed under this test category.
RapidTest	✓				
Consumable	✓				
DS ABENGOUROU-> CAT ABENGOUROU					
CD4	✓	✗	✗	0	Site has not instrument listed under this test category.
Chemistry	✓	✗	✗	0	Site has not instrument listed under this test category.
Hematology	✓	✗	✗	0	Site has not instrument listed under this test category.
VirallLoad	✓	✗	✓	1	
OtherTest	✓	✗	✗	0	Site has not instrument listed under this test category.
RapidTest	✓				
Consumable	✓				
DS ABENGOUROU-> CHR ABENGOUROU					
CD4	✓	✗	✓	1	
Chemistry	✓	✗	✓	2	
Hematology	✓	✗	✓	1	
VirallLoad	✓	✗	✓	1	
OtherTest	✓	✗	✗	0	Site has not instrument listed under this test category.
RapidTest	✓				
Consumable	✓				
DS ABENGOUROU-> CSU AFFALIKRO					
CD4	✓	✗	✗	0	Site has not instrument listed under this test category.
Chemistry	✓	✗	✗	0	Site has not instrument listed under this test category.
Hematology	✓	✗	✗	0	Site has not instrument listed under this test category.

Figure 105. Service and Platform tab

- 2 Double click a test for the site you want to edit.

The Site data entry form for the selected site is displayed.

- 3 See Editing Site Data on page 47 for more information on completing this form.

- 4 Click Save and Close to return to the Review your Data form.

Service and Platform information for the site will be updated on this form based on changes made.

### SUPPLY LIST

Supply List tab helps to identify products pack cost and usage rate.

- 1 Click Supply List tab to view the list of supplies included in the quantification.

The Supply List tab is displayed. (Figure 106)

Patient Behavior Service and Platform Supply List							
Product Name	Pack Size	Unit	Pack C...	Pac...	Usage ...	Quantify According To	Applied ...
<b>CD4</b>							
BD FACS Count-Control Kit	25	test	180.0000	1	1	CD4 Tests-FACSCount	Testing
BD FACS Count-FacsClean	5	litre	37.5000	1	1	CD4 Tests-FACSCount	Testing
BD FACS Count-FacsFlow	20	litre	32.5000	1	1	CD4 Tests-FACSCount	Testing
BD FACS Count-FacsRinse	5	litre	32.5000	1	1	CD4 Tests-FACSCount	Testing
BD FACS Count-Thermal Paper	1	roll	8.0000	1	1	CD4 Tests-FACSCount	Testing
FACS Count CD3/4 Reagent Kit	50	test	213.0000	1	1	CD4 Tests-FACSCount	Testing
<b>Chemistry</b>							
Fully (biosystem) Aiguille	1	Aiguille	153.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Calibrator 5x5m	25	ml	46.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Cleaner Bottle	500	ml	30.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Control Serum	25	ml	46.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Control Serum	25	ml	46.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Creatinine 4x5l	200	ml	55.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Cuvettes, 8x 50	3200	cuvette	55.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Glucose 1x200	200	ml	55.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Lampe	1	lampe	38.3000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Papier imprima	1	rouleau	15.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) Sample cups 1C	1000	cup	30.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) TGO/ASAT 1x2	200	ml	40.0000	1	1	Total patient samples run on Fi	Testing
Fully (biosystem) TGP/ALT 1x20l	200	ml	40.0000	1	1	Total patient samples run on Fi	Testing
<b>Hematology</b>							
Sysmex Kx 21 - Cell clean CL50	50	ml	67.8100	1	1	Hematology Tests-Sysmex KX	Testing
Sysmex Kx 21 - CellPack PK	20 20	Litre	50.9900	1	1	Hematology Tests-Sysmex KX	Testing
Sysmex Kx 21 - Control Eight-TF	1	Tube	22.1000	1	1	Hematology Tests-Sysmex KX	Testing
Sysmex Kx 21 - Contrôle Eight-T	1	Tube	22.1000	1	1	Hematology Tests-Sysmex KX	Testing
Sysmex Kx 21 - Contrôle Eight-T	1	Tube	22.1000	1	1	Hematology Tests-Sysmex KX	Testing

Figure 106. Supply List tab

Once the data has been reviewed—

- 1 Click Calculate button to move to run the forecast. See Calculate Results on page 118.

## Editing Forecasts

With the Morbidity Forecast Methodology screen displayed—

- 1 Double click on the forecast in the Morbidity Forecast Methodology window that you want to edit.
- 2 Click on the Previous button, Next button and/or Go Home button to move to each of the Working Area pages that needs to be edited.
- 3 Enter the new data.
- 4 Click on the Next button to move to the Calculating Results screen.
- 5 See Calculating Results for information on calculating the forecast results.

## Deleting Forecasts

With the Morbidity Forecast Methodology screen displayed—

- 1 Select the forecast that you want to delete.
- 2 Click on the Delete button.

## Calculating Results

With the Morbidity Forecast Methodology screen displayed—

- 1 Click on the Calculate button once you review your data to let the model automatically calculate how many lab supplies are required for the entire quantification.
- 2 On the pop up dialog click on yes to confirm.
- 3 Click on ok once the forecast is completed.

When the calculation is finished, different reports and summary charts will be displayed.



Figure 107. Forecast Results

## Patient Assumptions

This report displays the number of patients on treatment and in pre-ART for each month of the forecast period.

## CD4 Assumptions

This report lists the key CD4 assumptions and displays the number of CD4 tests expected during the forecast period, broken up into existing patients in treatment (EP-T), existing patients in pre-ART (EP-Pre ART), new patients in treatment (NPT), new patients in pre-ART (NP Pre ART), symptom-directed tests (SDT), invalid tests and wastage (wastage), and buffer stock reagents used to run controls (BSRTC).

## Hematology, Viral Load, Chemistry, and Other Test Assumptions

This report displays the number of tests expected in the forecast period; broken up



- ## TO NAVIGATE BETWEEN THE TWO CATEGORIES –

- ## TO PRINT THE REPORTS –

- Each of the chart reports in the tool has a feature to Copy, Save as Image, Print and Zoom.

- TO VIEW THE PREVIOUS REPORTS –**

- ## Adjust Morbidity Supply Forecast

**1** Click on Adjust supply Procurement button.

[illegible]

**2** Double click on the Quantify in Stock column.

- 3** Type in Quantity in Stock for each product.
- 4** Click on the Save button to save the data, and return to the Dashboard screen.

Or, click on the Cancel button to disregard all changes and return to the Dashboard screen.







HT95FC

16/MAR

M

HT95FC

16/MAR

HT95FC

16/MAR

HT95FC

16/MAR

F

HT95FC